



**PtoleMedic System**

**USA-Instructions for Use**

**MRI Protocol Quick Reference Guide**

**Document ID: US-LB-72-01-005(A)**

**December 21, 2020**

# Lento Medical Innovation, Inc. Imaging Quick Reference Guide for 1.5T and 3T Scanners

These PtoleMedic System instructions for use are a quick reference guide for an MRI Technologist while setting up the listed MRI scanning equipment. The most current MRI set-up guide is always available on-line for reference or download. Lento Encourages technologists to verify the use of the most recent version frequently.



**WARNING:** Please note that the use of MRI scanners in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted in these patients.



**WARNING:** The use of the PtoleMedic System software in pediatric patients has not been studied, and the results of such use in these patients are unknown. We do not recommend that such surgery not be attempted.



**Manufacturer:**

Lento Medical Innovation, Inc.  
15110 Northwest Freeway, Suite 150  
Houston, TX 77040 USA  
+1 (510) 413-3230

# CE 2460

**The EU Authorized Representative for the PtoleMedic System:**



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany  
Tel.: +49-511-6262 8630  
Web: [www.mdss.com](http://www.mdss.com)

**Doc. #: US-LB-72-01-005 (A)**

**Issue Date: 12/21/2020**

**Page 2 of 4**

**CONFIDENTIAL**

This document is the property of Lento Medical Innovation, Inc. and shall not be copied or used without written uthorization.

| <b>3T MRI Scanners</b>   | <b>General Electric</b>  | <b>Philips</b>   | <b>Siemens</b>  | <b>Toshiba</b>  |
|--|--|--|---|---|
| <b>3 plane Localizer Scan Parameters:</b>                                | <b>GP FLEX (GE Users)</b><br>4mm x 1mm Skip,<br>24cm FOV,<br>Matrix 256 x192,                              | 4mm x 1mm Gap,<br>240mm FOV,<br>Voxel Size to equal<br>256 x192                                      | 4mm x 25% Distance<br>Factor,<br>240mm FOV,<br>Base Resolution<br>256 x 80% Phase<br>Resolution | 4mm x 1mm Spacing<br>24cm FOV,<br>Matrix 256x192  |
| <b>Coronal Knee:<br/>Pulse Sequence</b>                                  | FRFSE-XL<br>CORONAL  | TSE CORONAL  | TSF CORONAL   | TSF CORONAL   |
| <b>Mode</b>  | 2D   | 2D   | 2D  | 2D  |
| <b>Imaging Options</b>   | No Phase Wrap<br>ON,<br>Tailored RF,<br>3DGR (3D<br>Geometry<br>Correction for<br>software version<br>23+) | Fold Over<br>Suppression ON,<br>100% Sampling,<br>"Default" Selected<br>for Distortion<br>Correction | 100% Phase<br>Oversampling,<br>3D Distortion<br>Correction Filter ON                            | Fold Over<br>Suppression<br>ON,<br>IDC (Intelligent<br>Distortion Correction)<br>Selected |
| <b>TE (Echo Time)</b>  | Min Full   | ~24 to 35 (28<br>nominal)  | ~24 to 35 (28<br>nominal)   | ~24 to 35 (28<br>nominal)   |
| <b>TR (Repetition Time)</b>  | Use TR to get<br>series in one<br>acquisition  | Use TR to get the<br>shortest scan time  | Use TR to get the<br>shortest scan time   | Use TR to get the<br>shortest scan time   |
| <b>Flip Angle (Deg)</b>  | 90   | 90   | 120   | 90<br>Flop Angle 160  |
| <b>Echo Train Length (ETL),<br/>Turbo Spin Factor (TSF)</b>              | 7  | 8  | 7   | 7   |
| <b>FOV</b>   | 16cm   | 160mm  | 160mm   | 16cm  |
| <b>Slice Thickness (mm)</b>  | 3mm  | 3mm  | 3mm   | 3mm   |
| <b>Spacing/Skip/Gap/Distance<br/>Factor (mm)</b>                         | 0  | 0  | 0%  | 0   |
| <b>Scan Matrix/Voxel Size<br/>Base Resolution x Phase<br/>Resolution</b> | 256 x 160  | Voxel to equal<br>256 x 160  | 256 x 75% Base<br>Resolution  | 256 x 160   |
| <b>NEX/NSA/Averages</b>  | 2  | 2  | 1   | 2   |
| <b>Frequency Direction</b>   | S/I  | S/I  | H/F   | S/I   |

| <b>3T MRI Scanners</b>   | <b>General Electric</b>   | <b>Philips</b>   | <b>Siemens</b>  |
|--|---|--|---|
| <b>3 plane Localizer Scan Parameters:</b>                                | <b>GP FLEX (GE Users)</b><br>4mm x 1mm Skip, 24cm FOV,<br>Matrix 256 x192,  | 4mm x 1mm Gap,<br>240mm FOV,<br>Voxel Size to equal<br>256 x192                                    | 4mm x 25% Distance<br>Factor,<br>240mm FOV,<br>Base Resolution<br>256 x 80% Phase<br>Resolution |
| <b>Coronal Knee:<br/>Pulse Sequence</b>                                  | FRFSE-XL CORONAL  | TSE CORONAL  | TSF CORONAL   |
| <b>Mode</b>  | 2D  | 2D   | 2D  |
| <b>Imaging Options</b>   | No Phase Wrap ON,<br>TRF (Tailored Radio<br>Frequency),<br>3DGR (3D Geometry<br>Correction for software<br>version 23+) | Fold Over<br>Suppression R/L,<br>100% Sampling,<br>"Default" Selected for<br>Distortion Correction | 100% Phase<br>Oversampling,<br>3D Distortion<br>Correction Filter ON                            |
| <b>TE (Echo Time)</b>  | Min Full  | ~24 to 35 (28<br>nominal)  | ~24 to 35 (28<br>nominal)   |
| <b>TR (Repetition Time)</b>  | Use TR to get series in<br>one acquisition or<br>shortest scan time   | Use TR to get the<br>shortest scan time  | Use TR to get the<br>shortest scan time   |
| <b>Flip Angle (Deg)</b>  | 90  | 90   | 120   |
| <b>Echo Train Length (ETL),<br/>Turbo Spin Factor (TSF)</b>              | 7   | 8  | 7   |
| <b>FOV</b>   | 16cm  | 160mm  | 160mm   |
| <b>Slice Thickness (mm)</b>  | 3mm   | 3mm  | 3mm   |
| <b>Spacing/Skip/Gap/Distance<br/>Factor (mm)</b>                         | 0   | 0  | 0%  |
| <b>Scan Matrix/Voxel Size<br/>Base Resolution x Phase<br/>Resolution</b> | 256 x 160   | Voxel to equal<br>256 x 160  | 256 x 75%Base<br>Resolution   |
| <b>NEX/NSA/Averages</b>  | 2   | 2  | 1   |
| <b>Frequency Direction</b>   | S/I   | S/I  | H/F   |