



## PtoleMedic System

# Instructions for Use – MRI Protocol Quick Reference Guide

Document ID: LB-72-01-005

Model: Revision C

April 12, 2022

*Author:*

Name: **Kenneth Jung**

**4/12/2022**

Date

*Reviewed and Approved by:*

Name: **David Schlerf**

**4/12/2022**

Date

Doc. #: LB-72-01-005 (C)

Issue Date: 4/12/2022

Page 1 of 4

**CONFIDENTIAL**

This document is the property of Lento Medical Innovation, Inc. and shall not be copied or used without written authorization.

# Lento Medical Innovation, Inc. Imaging Quick Reference Guide for 1.5T and 3T Scanners

This is a quick reference guide for the MRI Technologist as an aid in rapidly setting up the listed MRI scanning equipment. There is a full MRI reference guide for use in positioning and patient orientation which is available upon request if it has not already been supplied or is misplaced.



**WARNING:** Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.



**WARNING:** The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery not be attempted.



**Manufacturer:**

Lento Medical Innovation, Inc.  
15110 Northwest Freeway, Suite 150  
Houston, TX 77040 USA  
+1 (510) 413-3230

The following organization is the Authorized Representative for the PtoleMedic System:

# C € 2460



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany

<b>3T MRI Scanners</b>	<b>General Electric</b>	<b>Philips</b>	<b>Siemens</b>	<b>Toshiba</b>
<b>3 plane Localizer Scan Parameters:</b>	<b>GP FLEX (GE Users)</b> 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution	4mm x 1mm Spacing 24cm FOV, Matrix 256x192
<b>Coronal Knee: Pulse Sequence</b>	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL	TSF CORONAL
<b>Mode</b>	2D	2D	2D	2D
<b>Imaging Options</b>	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON	Fold Over Suppression ON, IDC (Intelligent Distortion Correction) Selected
<b>TE (Echo Time)</b>	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
<b>TR (Repetition Time)</b>	Use TR to get series in one acquisition	Use TR to get shortest scan time	Use TR to get shortest scan time	Use TR to get shortest scan time
<b>Flip Angle (Deg)</b>	90	90	120	90 Flop Angle 160
<b>Echo Train Length (ETL), Turbo Spin Factor (TSF)</b>	7	8	7	7
<b>FOV</b>	18cm	180mm	180mm	18cm
<b>Slice Thickness (mm)</b>	3mm	3mm	3mm	3mm
<b>Spacing/Skip/Gap/Distance Factor (mm)</b>	0	0	0%	0
<b>Scan Matrix/Voxel Size Base Resolution x Phase Resolution</b>	256 x 256	256 x 256	256 x 256	256 x 256
<b>NEX/NSA/Averages</b>	2	2	2	2
<b>Frequency Direction</b>	S/I	S/I	H/F	S/I

<b>3T MRI Scanners</b>	<b>General Electric</b>	<b>Philips</b>	<b>Siemens</b>
<b>3 plane Localizer Scan Parameters:</b>	<b>GP FLEX (GE Users)</b> 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
<b>Coronal Knee: Pulse Sequence</b>	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL
<b>Mode</b>	2D	2D	2D
<b>Imaging Options</b>	No Phase Wrap ON, TRF (Tailored Radio Frequency), 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression R/L, 100% Sampling, "Default" Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON
<b>TE (Echo Time)</b>	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
<b>TR (Repetition Time)</b>	Use TR to get series in one acquisition or shortest scan time	Use TR to get shortest scan time	Use TR to get shortest scan time
<b>Flip Angle (Deg)</b>	90	90	120
<b>Echo Train Length (ETL), Turbo Spin Factor (TSF)</b>	7	8	7
<b>FOV</b>	16cm	160mm	160mm
<b>Slice Thickness (mm)</b>	3mm	3mm	3mm
<b>Spacing/Skip/Gap/Distance Factor (mm)</b>	0	0	0%
<b>Scan Matrix/Voxel Size Base Resolution x Phase Resolution</b>	256 x 256	256 x 256	256 x 256
<b>NEX/NSA/Averages</b>	2	2	2
<b>Frequency Direction</b>	S/I	S/I	H/F