## **PtoleMedic System**

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COMPLETE

Patient

First Name	Middle Name	Last Name	Surgery Type	Surgery Side	Default Alignment	Date of Purchase Order
Friday	F	Sunday	TKR	Left	Neutral Boundary	Jul 27, 2020
Date of Birth	Sex	Age	Implant Information	Ship To	Bill To	
Apr 05, 1958	male	62	SHITHANEPIH GENESISII	July Distributor	July Distributor	

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
User:	Aug, July
Office Name:	July MRI Center
Office Address:	Third Street Berlin, 7454654, DE
Email Address:	won4@lento-medical.com
Phone Number:	49132132132132
Fax Number:	4

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Surgeon List							
Last Name	First Name	Middle Name	Office Name	Office Address	Email Address	Phone	Fax
Dr. Tuesday	Monday	H	Monday Office	1254 Humble [ Berlin, DE ]	won1@lento-medical.com	4902495000	1

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Case List							Search
Case #	Surgeon Name	Patient Name	DOPO	DOS	Surgery Type	Surgery Side	Status
K41en6L	Tuesday, Monday	Sunday, Friday	Jul 27, 2020	Sep 25, 2020	TKR	Left	Ready to Upload MRI

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MRI Upload Ready

**Surgeon**

First Name  
 Last Name  
 Middle Name  
 Title  
 Email Address  
 Phone

**Patient**


First Name  
 Last Name  
 Middle Name  
 Title  
 Email Address  
 Phone  
 Date of Birth  
 Date of Surgery  
 Surgery Type

**Drop Zone**  
 Click and Drag Patient MRI Files Here

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Surgeon:	Seidel, Markus	<a href="#">NEW ORDER</a> <a href="#">CURRENT ORDERS</a> <a href="#">COMPLETED ORDERS</a> <a href="#">VIEW PREFERENCES</a>
Office Name:	Dr. Markus Seidel	
Office Address:	Kastanienallee 2 Berlin, 10435, DE u@lentomedical.com	
Email Address:	u@lentomedical.com	
Phone Number:	49304483568	
Fax Number:		

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Case #	Patient Name	Order Date	Surgery Date	Surgery Type	Surgery Side	Status
K49g2R	Park, Se Jung	Jul 27, 2020	Nov 17, 2020	TKR	Right	Surgeon Approval Ready

Case ID: K49g2R  
Status: Surgeon Approval Ready

[COMPLETE](#)

Patient


First Name	Se Jung	Middle Name		Last Name	Park	Surgery Type	TKR	Surgery Side	Right	Default Alignment	Mechanical Alignment	Date of Purchase Order	Jul 27, 2020
Date of Birth	Nov 27, 1988	Sex	Female	Age	32	Item of Surgery	SHOULDER	Imped Information	One to	One to	One to	One to	One to

Surgeon	Distributor	Hospital	MRI Center
Name	Company Name	Hospital Name	Center Name
Seidel, Markus	Feder express	Charité universitätsmedizin	HOT-CONTRIA
Office Name	Contact Person	Contact Person	Contact Person
Dr. Markus Seidel	Express, Feder	Berlin, Hospital	Berlin, Berlin
Address	Address	Address	Address
Kastanienallee 2	Aden-von-Trost-Straße 3	Charitépl. 1	Karl-Liebk. Str. 23
Email Address	Email Address	Email Address	Email Address
u@lentomedical.com	u@lentomedical.com	u@lentomedical.com	u@lentomedical.com
Office Number	Company Number	Hospital Number	Center Number
49304483568	49304483568	49304483568	49304483568
Mobile Number	Mobile Number	Mobile Number	Mobile Number
49304483568	49304483568	49304483568	49304483568
Fax	Fax	Fax	Fax

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

#### CAUTION!

This is to confirm that you are rejecting the pre-planned adjusted measurements made on the surgeon approval page. Please fill out rejection form below.

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Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	10
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Ankle"

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
[1.5T 및 3T MRI 스캐너 Imaging reference]

**1. General Electric**

<b>3plane Localizer Scan Parameter</b>	GP FLEX (GE Users), 4mm x 1mm Skip 24cm FOV, Matrix 256x192
<b>Coronal Knee: Pulse Sequence</b>	FRFSE-XL CORONAL
<b>Mode</b>	2D
<b>Imaging Option</b>	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)
<b>TE (Echo Time)</b>	Min Full
<b>TR (Repetition Time)</b>	Use TR to get series in one acquisition
<b>Flip Angle (deg)</b>	90
<b>Echo Train Length (ETL)</b> <b>Turbo Spin Factor (TSF)</b>	7
<b>FOV (Field of View)</b>	18cm
<b>Slice Thickness(mm)</b>	2.5mm
<b>Spacing/Skip/Gap/Distance Factor(mm)</b>	0mm
<b>Scan Matrix/Voxel Size</b> (Base Resolution x Phase Resolution)	256 x 256
<b>NEX/NSA/Averages</b>	2
<b>Frequency Direction</b>	S/I

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\*Higher resolution than 256 x 256 is allowed, using software.

**2. Philips**

<b>3plane Localizer Scan Parameter</b>	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x 192
<b>Coronal Knee: Pulse Sequence</b>	TSE CORONAL
<b>Mode</b>	2D
<b>Imaging Option</b>	Fold Over, Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction
<b>TE (Echo Time)</b>	~24 to 35 (28 nominal)
<b>TR (Repetition Time)</b>	Use TR to get shortest scan time
<b>Flip Angle (deg)</b>	90
<b>Echo Train Length (ETL)</b> <b>Turbo Spin Factor (TSF)</b>	8
<b>FOV (Field of View)</b>	180mm
<b>Slice Thickness(mm)</b>	2.5mm
<b>Spacing/Skip/Gap/Distance Factor(mm)</b>	0mm
<b>Scan Matrix/Voxel Size</b> (Base Resolution x Phase Resolution)	256 x 256
<b>NEX/NSA/Averages</b>	2
<b>Frequency Direction</b>	S/I

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\*Higher resolution than 256 x 256 is allowed, using software.

### **3. Siemens**

<b>3plane Localizer Scan Parameter</b>	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
<b>Coronal Knee: Pulse Sequence</b>	TSE CORONAL
<b>Mode</b>	2D
<b>Imaging Option</b>	100% Phase Oversampling 3D Distortion Correction Filter ON
<b>TE (Echo Time)</b>	~24 to 35 (28 nominal)
<b>TR (Repetition Time)</b>	Use TR to get shortest scan time
<b>Flip Angle (deg)</b>	120
<b>Echo Train Length (ETL)</b> <b>Turbo Spin Factor (TSF)</b>	7
<b>FOV (Field of View)</b>	180mm
<b>Slice Thickness(mm)</b>	2.5mm
<b>Spacing/Skip/Gap/Distance Factor</b>	0%
<b>Scan Matrix/Voxel Size</b> (Base Resolution x Phase Resolution)	256 x 256
<b>NEX/NSA/Averages</b>	2
<b>Frequency Direction</b>	H/F

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Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	14
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Hip"

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