PtoleMedic System

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Case ID: K4fen6L
Status: Waiting for MRI

COMPLETE

Patient

First Name
Friday

Middle Name
F

Last Name
Sunday

Surgery Type
TKR

Surgery Side
Left

Default Alignment
Neutral Boundary

Date of Purchase Order
Jul 27, 2020

Date of Birth
Apr 05, 1958

Sex
male

Age
62

Date of Surgery
Sep 25, 2020

Implant Information
SHETHBHEPHEN GENESISII

Ship To
July Distributor

Bill To
July Distributor

Surgeon

Name
Tuesday, Monday

Office Name
Monday Office

Address
1254 Humble

Email Address
won1@lento-medical.com

Office Number
4902495800

Mobile Number
4902495800

Fax
1

Distributor

Company Name
July Distributor

Contact Person
Aug, July

Address
Houston, TX 77040, US

Email Address
won2@lento-medical.com

Company Number
11235464565

Mobile Number
11235464565

Fax
2

Hospital

Hospital Name
July Hospital

Contact Person
Aug, July

Address
Second Street

Email Address
won3@lento-medical.com

Hospital Number
49352435456456

Mobile Number
49352435456456

Fax
3

MRI Center

Center Name
July MRI Center

Contact Person
Aug, July

Address
Third Street

Email Address
won4@lento-medical.com

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
User:	Aug, July
Office Name:	July MRI Center
Office Address:	Third Street Berlin, 7454654, DE
Email Address:	won1@lento-medical.com
Phone Number:	49132132132132
Fax Number:	4

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Surgeon List							
Last Name	First Name	Middle Name	Office Name	Office Address	Email Address	Phone	Fax
Dr. Tuesday	Monday	H	Monday Office	1254 Humble [Berlin, DE]	won1@lento-medical.com	4952495000	1

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
Case List							Search
Case #	Surgeon Name	Patient Name	DOPO	DOS	Surgery Type	Surgery Side	Status
KallenB	Tuesday, Monday	Sunday, Friday	Jul 27, 2020	Sep 25, 2020	TKR	Left	Ready to Upload MRI

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Surgeon:	Seidel, Markus	NEW ORDER CURRENT ORDERS COMPLETED ORDERS VIEW PREFERENCES
Office Name:	Dr. Markus Seidel	
Office Address:	Kastanienallee 2 Berlin, 10435, DE u@lentomedical.com 49304483568	
Email Address:		
Phone Number:		
Fax Number:		

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Current Orders							Search
Case #	Patient Name	Order Date	Surgery Date	Surgery Type	Surgery Side	Status	
K49g2R	Park, Se Jung	Jul 27, 2020	Nov 17, 2020	TKR	Right	Surgeon Approval Ready	

Case ID: K49g2R
Status: Surgeon Approval Ready

COMPLETE

Patient


First Name	Se Jung	Middle Name		Last Name	Park	Surgery Type	TKR	Surgery Side	Right	Default Alignment	Mechanical Alignment	Date of Purchase Order	Jul 27, 2020
Date of Birth	Nov 27, 1989	Sex	Female	Age	40	Date of Surgery	Nov 17, 2020	Implant Information	STYKER TRIUMPH	Ship To	Charité Universitätsmedizin	Ship To	Fedex express

Surgeon	Distributor	Hospital	MRI Center
Name Seidel, Markus	Company Name Fedex express	Hospital Name Charité universitätsmedizin	Center Name HRT-Zentrum
Office Name Dr. Markus Seidel	Contact Person Friedrich, Fedex	Contact Person Berlin, Hospital	Contact Person Hri, Berlin
Address Kastanienallee 2	Address Aden-von-Troitz-Strasse 1	Address Charitépl. 1	Address Kallier Str. 23
Email Address u@lentomedical.com	Email Address u@lentomedical.com	Email Address u@lentomedical.com	Email Address u@lentomedical.com
Office Number 49304483568	Company Number 49304483568	Hospital Number 49304483568	Center Number 49304483568
Mobile Number 49304483568	Mobile Number 49304483568	Mobile Number 49304483568	Mobile Number 49304483568
Fax	Fax	Fax	Fax

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

CAUTION!

This is to confirm that you are rejecting the pre-planned adjusted measurements made on the surgeon approval page. Please fill out rejection form below.

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Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	10
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Ankle"

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[1.5T 및 3T MRI 스캐너 Imaging reference]

1. General Electric

3plane Localizer Scan Parameter	GP FLEX (GE Users), 4mm x 1mm Skip 24cm FOV, Matrix 256x192
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL
Mode	2D
Imaging Option	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)
TE (Echo Time)	Min Full
TR (Repetition Time)	Use TR to get series in one acquisition
Flip Angle (deg)	90
Echo Train Length (ETL) Turbo Spin Factor (TSF)	7
FOV (Field of View)	18cm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

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*Higher resolution than 256 x 256 is allowed, using software.

2. Philips

3plane Localizer Scan Parameter	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x 192
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	Fold Over, Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	Use TR to get shortest scan time
Flip Angle (deg)	90
Echo Train Length (ETL) Turbo Spin Factor (TSF)	8
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

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*Higher resolution than 256 x 256 is allowed, using software.

3. Siemens

3plane Localizer Scan Parameter	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	100% Phase Oversampling 3D Distortion Correction Filter ON
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	Use TR to get shortest scan time
Flip Angle (deg)	120
Echo Train Length (ETL) Turbo Spin Factor (TSF)	7
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor	0%
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	H/F

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Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	14
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Hip"

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