PtoleMedic System

Instructions for Use

Document ID: LB-72-01-006
Model: Revision H

October 2, 2026

Author:

Name: **Alex Park** **10/2/2026**
Date

Reviewed and Approved by:

Name: **Won Choi** **10/2/2026**
Date

TABLE OF CONTENTS

1. Instructions for Use - Adding a New Case…… 7
2. Instructions for Use - Uploading MRI to Database (MRI Center)……15
3. Instructions for Use - Template Design Approval (Surgeon)…….21
4. Instructions for Use – Clinical MRI Scanning……31
5. Instructions for Use – MRI Protocol Quick Reference Guide……48
6. Instructions for Use – T.K.R. Surgical Procedure……52
7. Instructions for Use - Hospital Cleaning, Sterilization, and Disposal……76

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Introduction

PtoleMedic System *web-based software* offers a simplified Total Knee Replacement (TKR) surgery planning interface. Surgical planning software is widely available but is generally PC or local server based. The use of such programs tends to be time consuming, expensive and requires frequent updating and/or IT maintenance. The PtoleMedic System offers **"log-on-and-use"** surgical planning tools and custom-made cutting guide creation. And the system is **low-to-no maintenance** requiring no IT support to install, use or maintain. The PtoleMedic® System web-based software provides direct, focused assistance on setting angles, rotation, cut plane, offsets and implant sizing as appropriate to TKR surgery.

Warning:

- The PtoleMedic System has not been evaluated in a pediatric population therefore performance in such cases is unknown and not recommended.
- The PtoleMedic System software is a prescription only medical device.
- The software is not a substitute for critical thinking and intra-operative adjustment of surgical goals based on education, training and experience of the practitioner.
- Surgical planning software only provides and documents useful alignment and orientation information based on specific individual anatomic data obtained from current MRI image sources.
- The software does not provide an absolute or only solution plan for joint replacement surgery, it only documents one possible approach, no surgical philosophy is recommended.
- The use of the software on mobile platforms, such as Google Android, Apple iPhone, or similar tablet-type devices has not been validated and is not recommended at any time.
- These are patient specific, single use, disposable instruments.
- Do not attempt to reuse, recondition or re-sterilize.
- Do not alter the instruments in any way.
- PtoleMedic System guides are to be used by a surgeon trained in the use of personalized instrument surgery (custom guides).
- Be aware that these patient specific instruments have been manufactured based of MRI scans of the patient. If the patient's anatomy has changed significantly since the time of the MRI scan, the instruments should not be used.
- The instruments should be properly cleaned before sterilization. Do not use if the instruments are broken, cracked, or when loose powder is present.
- The instruments in the package are provided non-sterile.
- The instruments in the package must be sterilized prior to use.
- Federal Law (USA) restricts this device to sale by or on the order of a physician.

Cautions:

- The use of aged (>3 months) MRI image files is not recommended. Accuracy of planning and guide fit will diminish with evolving or changing disease processes.

Precautions:

- Use only MRI data of recent origin obtained per established PtoleMedic System designated MRI protocols.
- Care should be taken to minimize excessive heat buildup due to friction between the PSI instruments and other instrumentation, such as drills. Excessive heat buildup can lead to deformation of the PSI instrument.
- Do not place heavy instruments on top of the cutting guides instruments during sterilization.

Limitations

- Metallic implants in or near the affected joint are known to interfere with the MRI images.
- Provides estimate of implant sizing only. Optimal implant size can only be determined during surgery and may differ from sizes projected during planning. Most estimated implant sizes will typically fall within one size of estimation.
- The software is not for use in planning revision/replacement surgery in persons already having implants in the affected joints.
- DICOM to JPEG file creation – The use of planning images based on the JPEG image format limits some file attributes encoded with the original DICOM images. For instance, some information contained in the DICOM header including image manipulation functions like pan, zoom, full gray scale control, contrast adjustment and pixel resizing are limited or unavailable.
- Digital x-ray data is not acceptable for planning purposes; the files must be MRI images.

Contra-indications

- Do not use in patients with active infection of the knee joint area.

User System Requirements

1. High speed internet access
2. A color monitor or display with VGA image display capabilities with 1280 x 800 resolution (recommended)
3. A mouse, touch pad or similar input device to interact with the surgery case
4. A keyboard (virtual or physical) for text input
5. A PC/Workstation w/1GB of RAM or more capable of running Windows or Mac OS
6. A Web Browser, (Internet Explorer, Firefox, Google Chrome or Safari)

Clinical Benefit

PtoleMedic System enhances the accuracy of implant placement, potentially leading to better outcomes in surgery and improvement in quality of patient life.

Indication for Use
The PtoleMedic System is an on-line orthopedic surgical planning software system. MRI images supply data sufficient to allow accurate modeling of anatomy for on-line surgical planning before knee joint replacement surgery. The surgeon preoperatively plans, reviews, adjusts and orients implant images to select implant size and create an idealized surgical plan for the first bony cuts only.

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively

Intended Users:
These instructions are intended to guide Surgeons, Physicians, and MRI Technologist in using the PtoleMedic System software for creating new patients, MRI scanning of patients, Uploading patient MRI scans, and approving template designs.

NOTICE: ***Software Anomalies***. The PtoleMedic System software supports the following Web Browsers, Internet Explorer (Ver 9.0) or newer, Firefox (Ver 15.0) or newer, Google Chrome (Ver 25.0) and Safari (Ver 4.0) or newer. Such browsers may be used with various operating systems, including those supporting Windows OS (XP, 7, Vista) or newer or Mac OS X (Ver 10.0) or newer.
Depending on your browser or operating system your screen may look slightly different depending on the software environment. This may include available fonts, spacing between graphical user interface (GUI) components, and some other minor appearance items.

NOTICE: The PtoleMedic System Web software does not interact directly with DICOM generating equipment (PACS devices). The image files needed by the PtoleMedic System Web software originate from the MRI/CT center as DICOM files which are transferred to the company server by Imaging Center personnel using DICOM compliant processes. File transfer protocols or file transfer requests are not generated or provided.

NOTICE: The PtoleMedic System Surgical Planning Software program does not interact with or affect DICOM image source equipment (PACS devices).

NOTICE: In the event of a serious incident in relation to the PtoleMedic System, report the serious incident to Lento Medical Innovation, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Updates: For updates to this document, please contact the Help Desk support team.
Help Desk: If this document does not address your question please contact Lento Medical Innovation, Inc.
Telephone: The Help Desk and Customer Service support teams can be reached via:
United States and Canada during normal operating hours +1 (510) 413-3230.

Company URL: http://www.lentomedical.com
PtoleMedic System URL: http://www.lentomedical.net

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that the use of MRI scanners in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be performed, please verify the MRI order with the physician.

WARNING: Making custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. For this reason, it is generally recommended that such patient not be entered into the patient database. Please verify this request with the ordering physician.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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PtoleMedic System

Instructions for Use - Adding a New Case

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively

Intended Users:
This is a reference guide for the office manager, surgical planner/scheduler or physician as an aid in smoothly registering a new patient into the PtoleMedic System database.

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that the use of MRI scanners in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be performed, please verify the MRI order with the physician.

WARNING: Making custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. For this reason, it is generally recommended that such patient not be entered into the patient database. Please verify this request with the ordering physician.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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

Introduction

These instructions are intended to guide the surgeon or surgeons' representative in creating a new case for a patient. Adding a new case through the Lento Medical Innovation web portal (www.lentomedical.net) will create a new order for the cutting guide.

1. **Log In**
 Lento Medical Admin or Sales personnel will assist the surgeon with registering an account for the Lento Medical Innovation web portal. Once a username and password are created, the surgeon will use the information to log in.
 i. Enter Username
 ii. Enter Password
 iii. Click "Login" button


2. **Surgeon Landing Page**
 The surgeon landing page will appear once logged in.

3. Adding a New Case

From the Surgeon Landing page, click the “NEW ORDER” button to create a new case.


As seen in the above image, select surgery type by clicking “TKR”.

If the Surgeon User already has more than one pre-set preference, then user can select a preference (maximum of 3 preferences).

4. Entering Information for New Case

A. Add New Patient

A new page will be presented to add the following information:


I. First Name*
II. Middle Name

III. Last Name*
IV. Sex*
V. Date of Birth*
VI. Date of Surgery*
VII. Surgery Side*


Note: * fields are mandatory

B. Check Preference Information

I. Check distributor information
II. Check Hospital information
III. Check MRI Center information
IV. If any of the information are incorrect, STOP contact your distributor or Lento admin.
V. If you want to change the preference, go back to the previous screen and select a different preference.
VI. If all the information is correct, move forward to the next step.

C. Confirm Surgery Preference Information

I. Billing
 i. Billing party pre-set, (Hospital, Distributor, or Surgeon)
II. Shipping
 i. Ship to location preset (Hospital, Distributor, or Surgeon)
III. Implant
 i. Implant Brand pre-set.
IV. Alignment Default Position
 i. Select alignment default position (Mechanical or Natural Boundary).
V. Click “Submit”, to submit new order for case creation.


If you want to change the preference, go back to the previous screen and select a different preference.

5. New Case Creation

Once a case is created, this will create a new order in the PtoleMedic System based on the information provided. A unique serial number, the “case ID number”, is assigned to the new patient case as seen below.

Sample Case:


If the new case created requires any changes or updates, please contact a Lento Medical personnel or the sales representative for further assistance.


I. Collapse button will hide preference information to minimize window space.


6. Tutorial Video

i. You can also access a tutorial video link located at the top right menu for a visual demonstration of how to create a new order.

HOME LOGOUT
Tutorial


Surgeon: Tuesday, Monday
Office Name: Monday Office
Office Address: 1254 Humble
Berlin,
10719, DE
Email Address: won1@lentomedical.com
Phone Number: 4902495800
Fax Number: 1
NEW ORDER
CURRENT ORDERS
COMPLETED ORDERS
VIEW PREFERENCES

Tutorial
Surgical Procedure
Surgeon Console

Distributor
Hospital
Surgeon Console
LENTO
MEDICAL INNOVATION
Lento Surgeon Console
Surgeon Training
Revision July 2020
0:01 / 2:35
Distributor
Distributor

7. Logout

If the surgeon user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.

PtoleMedic System

Instructions for Use - Uploading MRI to Database (MRI Center)

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:
This is the reference guide for the MRI Technologist as an aid in uploading the MRI scanning images to the Lento Medical Innovation cloud database. There is a separate full MRI reference guide for use in positioning and patient orientation. Also, a quick MRI setting reference guide is available for use by experienced technologists and both are available upon request if they have not already been supplied or is misplaced.

MD

UDI

Basic UDI-DI: ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery be attempted.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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
Introduction

These instructions are intended to guide the MRI Center to upload MRI files through the Lento Medical Innovation web portal (www.lentomedical.net).

1. **Log In**

 Lento Medical Admin or Sales personnel will assist the MRI representative with registering an account for the Lento Medical Innovation web portal. Once a username and password are created, the surgeon will use the information to log in.

 i. Enter Username
 ii. Enter Password
 iii. Click “Login” button


2. **MRI Center Landing Page**

 User landing page will appear listing all surgeons. See below.


 Click “Surgeon Orders” button to see a list of surgeons that have cases waiting for MRI files to be uploaded, as seen below. User will only see Surgeons in the list that have open orders waiting for MRI Files.


3. MRI Uploading Page

From the Surgeon list, User will the click the surgeon's last name, which will direct the User to Case List page that shows all cases that are waiting for MRI files to be uploaded for the respective surgeon. The status of the case will show as "Ready to Upload MRI", see below image. NOTE: Before uploading MRI files, verify the patient information and the Case ID (unique serial number):

User will click the case number and the MRI upload page will appear, see below.

Note: before uploading the MRI file, put all the MRI files "Coronal Kee", "Coronal Hip", and "Coronal Ankle" into one folder.

4. Uploading MRI

There are two ways you can upload MRI files.

A. Drag and drop
 User can drag the MRI file folder and drop the file folder into the "Drop Zone"

B. Manual folder selection
 User can click the "Drop Zone" to select the folder and click the "Upload" button.


5. **Uploading Progress**

 Once the User as begun uploading the MRI files, there will be an uploading progress bar as seen below.

6. **MRI upload confirmation**

 Once MRI files are uploaded correctly and successfully, a confirmation popup will present itself as seen below.


 Once you click “OK”, user will be re-directed to user landing page.

 Note: once MRI is uploaded successfully, case will no longer be on the case list.

 Repeat from Step 2 to upload MRI files for a different case.

7. **Logout**

 If the user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.



8. **Follow up**

 Our technican will review the MRI files that had being upload and contact you if there is any issues.

9. **Tutorial Video**

 i. You can also access a tutorial video link located at the top right menu for a visual demonstration of how to upload MRI Files.

PtoleMedic System

Instructions for Use - Template Design Approval (Surgeon)

Intended Purpose:

The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:

These instructions are intended to guide the surgeon in the surgical planning and template design, review and approval stage for an order of the PtoleMedic System guides. Approval is provided online through the Lento Medical Innovation web portal (www.lentomedical.net).

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC

UDI-DI/PI: +G283PMS0010/&&71E

SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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
1. **Log In**

 Lento Medical Admin or Sales personnel will assist the surgeon with registering an account for the Lento Medical Innovation web portal. Once a username and password is created, the surgeon will use the information to log in.

 i. Enter Username
 ii. Enter Password
 iii. Click “Login” button


2. **Surgeon Landing Page**

 The surgeon landing page will appear listing basic information provided by the Surgeon. To see active cases, press the “Current Orders” button.


 The status of a case that requires surgeon approval will show “Surgeon Approval Ready”


3. View Surgeon Planning Reference

The Surgeon can click the case ID number (unique serial number) and the patient case detail page will appear.

From the patient case detail page, surgeon can view case information.

Press “EXPNAD” button to view preference details shown below.


Press “COLLAPSE” to hide the preference.


Press the “Surgeon Approve” button and the Surgeon will be directed to the “Surgeon Template Design Approval” Page

4. Surgeon Template Design Approval Page

On this page, surgeon will have the option to view each coronal, sagittal, and axial view of the MRI surgical planning results.


Each Coronal, Axial, and Sagittal view will be enlarged by double clicking the corresponding MRI image, shown as below.


Each arrow button will move the implant positioning (left, right, up, and down)
The rotate icon/button will rotate the implant by 3 degrees.


The overall limb alignment may be adjusted by clicking the right or left arrow on either side of the knee graphic image. The circle at the top of the graphic represents the center of the femoral head and the brackets at the bottom represent the center of the ankle joint. Each right or left arrow click adjusts the varus/valgus angle of the knee 0.5 degree or 0.5mm, max +- 5 degree or mm. The knee image is moved medially or laterally while maintaining the upper and lower boundary area within the femoral head and the median of the ankle. The resection lines of the rectangular resection box remain parallel to the floor.


Hip button will display Femoral Head Center in pop up.

Femur button will display Distal Femur Center in pop up.

Tibia button will display Proximal Tibia Center in pop up.

Ankle button will display Ankle Center in pop up.


After the surgeon has finished reviewing the planning data, at the bottom there will be four action buttons displayed. A button labeled reset, a button labeled reject, button labeled accept, and a button labeled close.

If the user presses the close button, system shall discard any saved changes and restore the original values and close the surgeon template approval page. The page will be redirected to the case detail page.

If the user presses the reset button, all the values will revert back to the original values from the beginning prior to any adjustments made by the surgeon user.

If the user presses the reject button, the system shall reject the entire plan. The surgeon will need to submit a reason for rejecting the plan as shown below.


Press Ok to submit the reason for rejection.


If the user presses the accept button, the OMM shall accept the current planning values. The user shall not be allowed to change the values after acceptance.

5. Tutorial Video

i. You can also access a tutorial video link located at the top right menu for a visual demonstration on how to approve the Template Design.


6. Logout

 If the surgeon user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.

PtoleMedic System

Instructions for Use – Clinical MRI Scanning

This is a quick reference guide for the MRI Technologist as an aid in rapidly setting up the listed MRI scanning equipment. There is a full MRI reference guide for use in positioning and patient orientation which is available upon request if it has not already been supplied or is misplaced. There is also an MRI upload reference guide to assist with uploading completed images to the Lento Medical Innovation cloud database.

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery not be attempted.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer of Record:
Kalms Consutling GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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A. PATIENT SET-UP

1. Patient should be metal/implant free in the scan volume.
2. Place patient supine, feet first, for lower extremity with affected limb as close to IsoCenter as possible
3. Consider using small foam wedges to stabilize or support the joint.
4. Use tape, sponges and/or Velcro straps on extremity to maintain anatomic “toes up” position.
5. Remind patient to remain motionless during the entire scan.
6. Landmark or Zero table
7. Begin scanning at the Knee series and move to Distal Tibia (Ankle) for the Coronal T1 Series and complete the study with the Hip series (see Figure 11).
8. Phased array coils (receive only coils) or the “Body Coil” may be used to scan the hip and ankle without removing the knee surface coil. New GEM suite GE systems have built in phased array coils in table.


B. SCANNER SET-UP

1. Enter patient name (Last, Middle, First)
2. Enter patient DOB. (dd/mm/yyyy)
3. Enter scan date (dd/mm/yyyy)
4. Enter left or right joint (Right/Left)
5. Enter patient gender (M/F)
6. Enter name of Surgeon (Last, Middle, First)
7. Enter name of imaging center

C. CORONAL FSE PD KNEE MRI PROTOCOL

1. Perform one 3 planes (Axial, Sagittal, Coronal) locator series.
2. Refer to "PtoleKnee Imaging Reference for 1.5T and 3T Scanners section for the Coronal scan parameters

D. CORONAL FSE PD KNEE SLICE POSITIONING

(a) **(b)**

Figure 2
Slice Positioning through Posterior Reference and Distal Reference Lines

1. Using the best axial and coronal locator images that show in Figure 2
 i. both posterior lateral and medial condyles, place slice parallel through the posterior reference for correct slice angle
 ii. Both distal lateral and medial condyles, place slice parallel through the distal reference for correct slice angle.


Figure 3
Slice Coverage

2. Slice coverage must include all the Femur including the condyles, Tibia and Patella as shown in Figure 3.
3. Adjust T1 T2 setting to show optimal image quality

Figure 4
Centering FOV

4. Using the best sagittal locator image that approximately shows the knee center, adjust the FOV to be centered at the “point” at the bottom of the distal femur as shown above in Figure 4.



Figure 5
Axial axis of the Knee joint

Not Acceptable **Acceptable**

5. Find the true Axial plane or axis of the knee joint as shown in Figure 5 by the yellow dash line. In addition, locate the sagittal MRI scout scan slice representing approximately the mid tibia, and tilt about 30 degrees. The upper bound of femur coverage must include femur anterior cortex or the beginning of the femoral shaft. The lower bound of tibial coverage must include the tibial shaft as much as possible, shown in Figure 6.


Figure 6
Slice Coverage

6. **Enter File Series Name: “Coronal Knee.”**

A. General Electric

3plane Localizer Scan Parameter	GP FLEX (GE Users), 4mm x 1mm Skip 24cm FOV, Matrix 256x192
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL
Mode	2D
Imaging Option	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)
TE (Echo Time)	Min Full
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	18cm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

B. Philips

3plane Localizer Scan Parameter	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x 192
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	Fold Over, Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	8
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

C. Siemens

3plane Localizer Scan Parameter	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	100% Phase Oversampling 3D Distortion Correction Filter ON
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700**)
Flip Angle (deg)	120
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor	0%
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	H/F

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

D. Toshiba

3plane Localizer Scan Parameter	4mm x 1mm Spacing, 24cm FOV, Matrix 256 x 192
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	Fold Over Suppression ON, IDC (Intelligent Distortion Correction) Selected
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90, Flop Angle 160
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	18cm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256x256 is optimal and preferred.

E. CORONAL T1 DISTAL/ANKLE MRI PROTOCOL

1. Perform one 3 plane (Axial, Coronal, Sagittal) locator scan
2. Perform a Coronal (T1FSE) series using the minimum scan
 i. parameters:

Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	10
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Ankle"

Figure 1
Coronal Slice Placement

F. CORONAL T1 DISTAL TIBIA/ANKLE SLICE POSITIONING

View the **Sagittal** locator to place the slice Coronally through the long axis of the Tibia. *The Axial locator slices are not used for alignment.* Center the FOV approximately one inch above the Tibial Talar joint space. See Figure 1 for slice position and FOV centering.

G. CORONAL T1 HIP MRI PROTOCOL AND SLICE POSITIONING

1. Perform one 3 plane (Axial, Sagittal, Coronal) locator scan.
2. Perform a Coronal (T1 FSE) series using the minimum scan parameters total below:
3. Place Coronal slice plane parallel through long axis of the femoral neck when viewing the axial locator as shown in Figure 8. Center the 14 slices to cover the Femoral Head, Neck, and Greater Trochanter as shown in Figure 9.


Figure 10
Center Hip FOV

Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	14
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Hip"

I. QUALITY CHECK AND METAL ARTIFACT INSTRUCTIONS

a. Coverage: Include 1 or 2 slices of soft tissue on both medial and lateral side of bony knee anatomy.
b. Motion: Make sure sequential transition of the bone/cartilage outline does not jump or shift back and forth between slices
c. If patient presents with metal in the affected lower extremity, follow these instructions
 - Metal Hardware Knee – Do **NOT** Scan the patient. Refer them back to the surgeon
 - Metal Hardware Hip and/or Distal Tibia – Use MARS (Metal Artifact Reduction Sequence) series. Keep in mind the FOV, Slice Thickness, and Spacing/Gap/Skip/Distance Factor must be set to our parameters. Do not use FAT SAT at any time.
 - For those facilities that do not use a MARS series, increase the bandwidth and NEX/NAQ/NSA to help increase signal to noise ratio.

II. FILE FORMAT & OUTPUT

a. Save all image series in DICOM format
b. Upload the following DICOM files to Lento Medical Innovation's database (www.lentomedical.net)
 i. **Coronal Knee**
 ii. **Coronal Hip**
 iii. **Coronal Ankle**



Hip Images scanned third
Knee Images scanned first
GLOBAL COORD.
Ankle Images scanned second

Figure 11


Scanning Order

I. KNEE, ANKLE, AND HIP MRI SAMPLE IMAGES


a. **KNEE MRI – file name: Coronal Knee**


(a)

(b)

(c)


(d)

(e)


(f)


(g)

b. **ANKLE MRI – file name: Coronal Ankle**


(a)

(b)


(c)

(d)


c. **HIP MRI – file name: Coronal Hip**


(a)

(b)

(c)

(d)

PtoleMedic System

Instructions for Use – MRI Protocol Quick Reference Guide

Lento Medical Innovation, Inc.
Imaging Quick Reference Guide
for 1.5T and 3T Scanners

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:
This is a quick reference guide for the MRI Technologist as an aid in rapidly setting up the listed MRI scanning equipment. There is a full MRI reference guide for use in positioning and patient orientation which is available upon request if it has not already been supplied or is misplaced.

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery be attempted.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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3T MRI Scanners	General Electric	Philips	Siemens	Toshiba
3 plane Localizer Scan Parameters:	**GP FLEX (GE Users)** 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution	4mm x 1mm Spacing 24cm FOV, Matrix 256x192
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL	TSF CORONAL
Mode	2D	2D	2D	2D
Imaging Options	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON	Fold Over Suppression ON, IDC (Intelligent Distortion Correction) Selected
TE (Echo Time)	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
TR (Repetition Time	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**
Flip Angle (Deg)	90	90	120	90 Flop Angle 160
Echo Train Length (ETL), Turbo Spin Factor (TSF)	7	8	7	7
FOV	18cm	180mm	180mm	18cm
Slice Thickness (mm)	2.5mm	2.5mm	2.5mm	2.5mm
Spacing/Skip/Gap/Distance Factor (mm)	0	0	0%	0
Scan Matrix/Voxel Size Base Resolution x Phase Resolution	256 x 256	256 x 256	256 x 256	256 x 256
NEX/NSA/Averages	2	2	2	2
Frequency Direction	S/I	S/I	H/F	S/I

3T MRI Scanners	General Electric	Philips	Siemens
3 plane Localizer Scan Parameters:	**GP FLEX (GE Users)** 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL
Mode	2D	2D	2D
Imaging Options	No Phase Wrap ON, TRF (Tailored Radio Frequency), 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression R/L, 100% Sampling, "Default" Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON
TE (Echo Time)	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
TR (Repetition Time	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**
Flip Angle (Deg)	90	90	120
Echo Train Length (ETL), Turbo Spin Factor (TSF)	7	8	7
FOV	16cm	160mm	160mm
Slice Thickness (mm)	2.5mm	2.5mm	2.5mm
Spacing/Skip/Gap/Distance Factor (mm)	0	0	0%
Scan Matrix/Voxel Size Base Resolution x Phase Resolution	256 x 256	256 x 256	256 x 256
NEX/NSA/Averages	2	2	2
Frequency Direction	S/I	S/I	H/F

PtoleMedic System

Instructions for Use – T.K.R. Surgical Procedure for the PtoleKnee Surgical Guide

PtoleMedic System
Instructions for Use
T.K.R. Surgical Procedure

IMPORTANT NOTE FOR THE SURGEON: PLEASE READ THE ENTIRE I.F.U. BEFORE USING THE ***PtoleKnee Surgical GUIDES*** CLINICALLY. THE SURGEON SHOULD BE FAMILIAR WITH THE PERSONALIZED SURGICAL INSTRUMENT CONCEPT.

Warnings:

- The ***PtoleMedic System*** has not been evaluated in a pediatric population; performance in such cases is unknown and not recommended.
- The ***PtoleMedic System*** is a prescription-only medical device.
- The ***PtoleMedic System*** is not a substitute for critical thinking and intra-operative adjustment of surgical goals based on the surgeon's education, training, and experience.
- ***PtoleMedic System*** only provides and documents functional alignment and orientation information based on specific individual anatomic data obtained from current M.R.I. image sources.
- The ***PtoleMedic System*** does not provide an absolute or only solution plan for joint replacement surgery; it only documents one possible approach.
- ***PtoleKnee Surgical Guides*** are one-time-use, disposable instruments.
- Do not attempt to reuse, recondition, or re-sterilize.
- Do not alter the custom guides in any way.
- ***PtoleKnee Surgical Guides*** are for use by a surgeon experienced in using personalized surgical instruments (PSI-customized guides).
- ***PtoleKnee Surgical Guides*** are patient-specific instruments planned and made based on M.R.I. scans for each named patient. New images should be obtained if the patient's anatomy/disease process has changed significantly since obtaining the M.R.I. scan.
- Examine the process and sterilize the ***PtoleKnee Surgical Guides*** before use. Do not use the guide if chipped, broken, cracked, or debris is present.

- The guides and their packaging are ***non-sterile***.
- ***PtoleKnee Surgical Guides*** may not be reused. They are for ONE-TIME-USE only.

Cautions:

- The use of aged (>3 months) M.R.I. image files is not recommended. Accuracy of planning and guide fit may diminish with evolving or changing disease processes.

Precautions:

- Use only M.R.I. data of recent origin obtained per established ***PtoleMedic System*** designated M.R.I. protocols.
- Take care to minimize excessive heat buildup from friction between P.S.I. instruments and other instrumentation, such as drills/saws. Excessive heat buildup can lead to debris or deformation of the ***PtoleKnee Surgical Guides*.**
- Do not place heavy instruments on top of the ***PtoleKnee Surgical Guides*** during sterilization.
- Surgical guides are made for scheduled surgery dates only! Storage conditions only require the guides to remain in their original packaging until processed for surgery. The shelf life of the guides is related to the patient's disease progression. If the patient's surgery is delayed for more than 3 months, the surgeon should determine if new M.R.I. data should be obtained and new guides made.

Limitations

- Metallic implants in or near the affected joint interfere with the M.R.I. images and may yield unreliable or useless images.
- The ***PtoleMedic System*** planning provides an estimate of implant sizing only. Exact implant size can only be determined during surgery and may differ from sizes projected during planning. Most estimated implant sizes fall within one size of estimation.
- The ***PtoleMedic System*** is not used in planning revision/replacement surgery in persons already having implants in the affected joints.
- Digital x-ray or C.T. data is not acceptable for guide production; the files must be M.R.I. images.

Contra-indications

- The ***PtoleKnee Surgical Guides*** should not be used when:
- Active infections of the knee or knee joint are present, Hip-Knee-Ankle alignment deformity greater than 6° Varus or valgus, or the case requires "tibia-cut-first" surgical techniques, cases of uni-condylar replacement, and for cases of T.K.R. revision surgery.

NOTICE: ***PtoleKnee Surgical Guides*** are intended to assist in executing a designated joint replacement (T.K.R.) surgery and must be created using the ***PtoleMedic System*** web application software. The guides are not reusable or transferable to any other person or surgery type.

1. Help Desk: For questions, please contact;
2. Lento Medical Innovation, Inc.

Telephone: United States telephone support +1 (510) 413-3230
Internet: If you have access to the Internet, you may reach the Online Help support:

Web: https://www.lentomedical.com or https://www.lentomedical.net

Importer of Record:
Kalms Consutling GmbH

Rheinstr. 45-46
12161 Berlin, Germany

Manufacturer:

3. **Lento Medical Innovation, Inc.**

15110 Northwest, Freeway, Suite 150
Houston, TX 77040 - U.S.A.
+1 (510) 413-3230
Help services are available free of charge to all registered customers.

Indications for Use:
The ***PtoleKnee Surgical Guides*** are for clinical use as a template or guide for use in total knee replacement orthopedic surgery to assist the surgeon in selecting or positioning orthopedic implants and guiding the marking of tissue before cutting or pining for a specifically named patient.

CAUTION: United States Federal law restricts this device to sale by or on the order of a Physician

Acknowledgments:
Lento Medical Innovation, Inc. acknowledges the assistance of the following orthopedic surgeon for his expertise, guidance, and time in the development of the surgical protocol:

Benjamin Soo-Il Song, M.D., Diplomat American Board of Orthopedic Surgery

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4. Introduction

These instructions cover the use of personalized surgical instruments designed using the ***PtoleMedic System*** surgical planning software to make T.K.R. femur and tibial component alignment guides.

Important Notes:
PtoleKnee Surgical Guide **Inspection**: Upon delivery, verify the presence and accuracy of the engraved information on the cutting guides.

- Serial # of *the* guides are in the format K, with five numerals, & an L or R (e.g., K12345L or K12345R)
- Patient's 1st initial & up to 10 letters of the last name
- Patient's Date of Birth

Indications/Contraindications: Review this entire section of all indications, contraindications, warnings, and precautions before ordering ***PtoleKnee Surgical Guides.***

Start with the femur: Lento does not recommend a tibia-first approach.

Incision Length: M.I.S. incisions are unlikely to provide sufficient access for the use of ***PtoleKnee Surgical Guides***. A standard incision is essential to ensure accurate placement of the ***PtoleKnee Surgical Guides***.

5. Pre-operative Imaging Scan:

The initial step in the ***PtoleMedic System*** process is a quality M.R.I. scan of the arthritic knee and a scout M.R.I. scan of the hip and the ankle. Lento recommends a 1.5 Tesla or higher magnet obtained from a qualified M.R.I. imaging center. **All scans should be obtained at least 14 calendar days before the surgical procedure date** and sent electronically to the Lento secure database via the ***PtoleMedic System*** M.R.I. web interface.

Before submitting scans for the ***PtoleMedic System*** knee procedure, pre-qualification of the imaging center by Lento is required. Strict scanning procedures and quality control measures apply to ensure accurate imaging of the patient's knee. Patients with a pacemaker, defibrillator, and large thigh circumference not fitting within the knee or torso coil or the ability to remain motionless for the scan are not recommended for the ***PtoleMedic System*** procedure.

6. The *PtoleMedic System*:

The ***PtoleMedic System*** is a Web-based approach to orthopedic surgical planning that enables surgeons to carefully preplan joint replacement procedures (T.K.R.) and personalize each patient's surgical parameters. Lento's proprietary software provides implant alignment and placement information based on the individual's anatomic data from medical M.R.I. images. It also allows the surgeon to request individualized custom surgical cutting guides for personalized alignment and positioning of the implants during surgery.

7. Sterilization of *PtoleKnee Surgical Guides*:

The Femur and Tibia cutting guides are supplied clean but **not sterile**. Hospital processing and sterilization recommendations accompany each guide shipment to the hospital or surgical center.

8. Recommended Saw Blade:

The recommended saw blade thickness is 1.27 mm (0.050 inch), and the saw blade length of 110-120 mm.

9.

10. Distal and Proximal Cutting Guides:

The distal femoral and proximal tibial cutting guides help set the Varus/valgus, flexion/extension (femoral), posterior slope (tibial), and proximal/distal positions of the planned femoral and tibial components.

Surgical Procedure Overview:
No special or unique surgical approach is required. The guides are functionally useful with almost all standard (non-MIS) surgical incisions and approaches for the knee.

Femur

Slide the femur guide from top to bottom in rolling mode until it fits the anterior condyles, gently pushing from the top down to firmly engage the femoral lateral guide arms. Check from anterior to posterior to ensure the guide sits flush against the distal condyles.

Hold the guide firmly by hand and secure it with two additional proximal pins inserted from the top.
Next, insert the angled pin.
Then drill the anterior holes to place the rotation pins.

Make sure not to remove soft tissue or osteophytes before positioning the guide. Osteophytes may be removed in case they interfere with the fitting of the guide.

Tibia

Insert the tibial guide from the top, checking the five contact points on the tibial surface, and slide it until the anteromedial step engages the anteromedial cortex.

Hold the guide firmly from the top.
First, insert the top posterolateral pins (collared and treated at the end).
Second, insert the anterior angled pin.
Third, insert the second anterior pin (more medial).

While holding the tibia from the top, make tibia cut starting from the medial side to posterior and lateral then take out the postero lateral pin to finish the tibia cut.

#	Steps	Procedure Overview	PST Instruments
1	*Knee Exposure*	Expose joint using standard approach. Allow enough exposure to allow seating of guide	NONE
2	*Excise A.C.L.*	Excise the ACL	NONE
3	*Femur Preparation*	Place Femur Guide anteriorly on the distal femur and slightly push down. (make sure the contact points are in place)	Femur Guide
4	*Femur Guide Placement*	1. Drill the two condyle pins first 2. Drill the angled third pin 3. Drill the frontal pin holes to set rotation (anterior)	Femur Guide
5	*Resect Distal Femur*	1. Verify guide placement and proceed with the Distal Cut. (Safely hold the guide in place with one hand) 2. Remove femur pins and guide	Femur Guide
6	*Tibial Exposure*	1.Expose tibia & excise meniscus & soft tissue at guide contact points	NONE
7	*Tibial Preparation*	1. Fix Tibial Guide on anterior surface and tibial plateau	Tibial Guide
8	*Tibial Guide Preparation*	1. Install posterolateral top pin (collared and treated at the end) 2. Drill the two anterior pins (check for stability and alignment)	Tibial Guide
9	*Resect tibia*	1. make tibia cut starting from the medial side to posterior and lateral then take out the postero lateral pin to finish the tibia cut. (Safely hold the guide in place with one hand) 2. Remove guide & pins	Tibial Guide
10	*Standard T.K.R.*	Return to standard T.K.R. procedure. (Use 4:1 block and Angel Wing to adjust pin rotation and check for notching on the femur).	NONE


The Surgical Procedure:

Warning: If the Case ID markings do not match the patient and each other, **do not Use** the ***PtoleKnee Surgical Guides*** for the surgery. **Use Standard Instruments**. Notify Lento as soon as possible after the completion of the case.

Incision:

A median parapatellar incision beginning one or two finger-widths superior to the patella and ending at the approximate location of the medial edge of the tibial tubercle is usually sufficient (see illustration below). However, depending on surgeon preference, the median parapatellar or variations such as mid-vastus, sub-vastus or lateral can be used as long as sufficient exposure results allow proper guide insertion. The smaller actual MIS-type incisions do not provide sufficient surgical exposure for guide insertion.

Parapatellar Incision

The initial step proceeds as typical for T.K.R. procedures, opening, dissecting, and removing adipose, Anterior Cruciate ligament, and capsule tissue. Patellar dislocation is done as you usually would.

Exposure:
With the knee moderately flexed, the medial synovium is released from the mid-point of the patella proximally to a point superior to the trochlear groove. The posterior patellar tendon fat pad is excised from the joint line to the tibial tubercle. With the knee flexed to approximately, 70 degrees retract the quadriceps muscle in the usual fashion to expose-to-expose the anterior femoral cortex. Displace the patella laterally as usual to obtain full exposure.

Osteophytes:
Removal of osteophytes is usually optional. ***PtoleKnee Surgical Guides*** are planned to affix in areas devoid of osteophytes. Their presence is generally not a detriment to the proper fit of the guides. If an osteophyte compromises the fitting of the guide, remove it and proceed with fitting the guide, ensuring that all contact points are properly in place.

Placing the Guides:
Place the guides sequentially, beginning with the femur and finishing with the tibia. Whenever guides are placed and before any bone cuts are made, it is essential for the surgeon to visually verify and or measure the proposed orientation and angle of the impending cut. Templates and guides are not a substitute for sound clinical judgment. If any potential cut is believed to be inappropriate, the use of the guide should be abandoned, and the standard manual instruments used without delay completing the case.


Examine the Femur Guide:
First, ascertain that you have received the correct guide. The face of the ***PtoleKnee Surgical Guide*** contains the case identification information (red oval). Engraved here are the Guide Serial Number, Patient initial, and ten alpha characters of the last name, D.O.B., and knee descriptor.

A.P. View – Face


Next, examine the physical properties of the guide.

P.A. View - Interior

The four primary guide contact points are identified in the picture above. The locations of the pinning holes are marked with black arrows. Scrutinize all features; they should be smooth and free of debris.

Medial - Lateral View



The lateral view above shows the cut slot indicated with a black arrow. Also depicted are the guide contact points. Also, note the contact points from the previous A.P. guide view. Guide contact point (2) is directly behind contact point (1). As explained previously, scrutinize all surfaces and features with particular attention to the cut slot; all should be smooth and free of debris.

Placing the Femur Guide:


Initial placement of the femur guide is done by orienting the guide and the condyle contact feet slightly above the mid-point of the condylar curve to allow the superior feet to make initial contact slightly distal of the epicondylar area (blue arrow). The face of the guide (red line) should approximate a twenty-degree angle relative to the long axis of the femur (black line).

The ***PtoleKnee Surgical Guide*** should rest initially at a modest angle superiorly as shown. The guide body is centered over and within the trochlear groove with distal condylar feet resting lightly on the distal condyle surfaces. Apply slight downward pressure to allow the guide to rotate distally and posteriorly; the guide can be felt to stick slightly as it locates the intended fixation point.


Once the guide is rotated posteriorly, the anterior and posterior feet of the guide make light contact with the anterior lateral and medial sides of the condyles and the distal condylar surfaces, as shown (blue arrows). When properly positioned, the guide's face (red line) indicates the intended orientation of the cut plain. The black line is perpendicular to the red line and typically approximates the long axis of the femur. The blue line indicates the intended cut plain and indicates the depth of the intended condylar resection.

From the medial side, verify the position of medical guide contact points (red circles). Each contact point should lightly touch the surface of the tissue without using excessive force or excessive direct pressure. A narrow gap (less than 1.0 mm) may sometimes occur along the guide contact point surface.



From the lateral side, in the same fashion as the medial side, verify lateral contact points. Both contact points should lightly touch the surface of the tissue without using excessive force or direct pressure. A narrow (less than 1.0 mm) gap may sometimes occur somewhere along the contact point surface.

There are four tissue contact points (2 medial and 2 laterals) and one visual indicator (5th point). Verify the correct placement location of the guide before pinning or making any cuts.

The fifth point (visual indicator) is located over the femoral notch and must be viewed as shown in the two following images. This fifth guide point is a machined notch in the form of a "V" shaped slot aligned so that a pin or rod placed within it points to the center of the femoral


head. When the guide is appropriately placed, the V notch is located over the mid-line of the condylar notch (center of the knee).


Verification of the potential Varus/valgus angle is made using a fixation pin placed within the "V" notch on the anterior surface of the guide, as shown. The pin should point to the center of the femoral head (black arrow) if the guide placement matches the surgical plan. After establishing the projected Varus/valgus alignment is as desired, the projected condylar cuts are assessed and may be completed.

Verification of the projected initial cut plane for the medial and lateral condylar cuts (red line) of the distal femur is performed. This cut plane is typically perpendicular to the long axis of the femur (black line) or as dictated by the surgical plan if flex/extension adjustments are planned.

Use a resection checker (angel-wing) to assess the angle and depth of each condylar cut to verify they are as expected. Angle verification is essential; as the initial cut sets the final Varus/valgus angle, flex/extension, and IR/ER of the femur implant. The thickness of each condyle cut should match the surgical plan values. If satisfactory, fixation of the guide can follow.

An additional guide placement verification is performed to assess that the initial cut plane will not ultimately result in a resection that leads to the implant's notch or gapping following the 4:1 (chamfer) cutting block.

The vertical black line above represents the intended cut plane for the distal femur. The horizontal red lines approximate the final anterior and posterior cut planes resulting from a 4:1 chamfer-cutting block. The anterior and posterior cut angles may be nearly parallel or divergent depending on implant brand, as shown. In either case, the final anterior cut plane of the 4:1 cutting block should skim the anterior surface of the femur (blue arrow).

Drill and set the Pins


Pre-drilling pilot holes for the pins obtain the most accurate position of the guides. It is recommended that the distal condylar pins should be drilled and set first. The holes are not drilled to maximum depth, only deep enough to pierce the cortical bone. Insertion of the pins after pre-drilling reduces the tendency of threaded or smooth pins to "wander or walk" before they bite in. For the femoral guide, begin by placing the 2 top pins most proximal to the femoral head, followed by the single angled pin. Then drill the anterior pins for rotation and AP femur setting; if additional stability is required, place one fixed anterior pin that will be removed last at the distal femur condyle. It is not recommended to start from the anterior femur pin, as the elbow can shift up or down the femoral guide.


With the fingers or thumb of the non-drilling hand, stabilize the guide by moderate pressure over the trochlear groove area to prevent unwanted movement, maintain pressure until both condylar pin holes are drilled.

1) Using an orthopedic **3.2 mm** drill, insert the bit into the guide hole before turning on the drill motor, then **drill the 2 top pins most proximal to the femoral head (white arrow),** taking care to minimize contact with the walls of the guide hole. Check the ***PtoleKnee Surgical Guide*** contact points to ensure the guide has not shifted, correct as necessary. Replace this drill with a **3.2 mm** bayonet or trocar point pin and impact or screw into place. Remove the lateral drill bit and replace it with a pin in like fashion.


2) **Proceed to drill the single oblique (angled) pin. (blue arrow)**. Check the ***PtoleKnee Surgical Guide*** contact points to ensure the guide has not shifted, correct as necessary. Replace this drill with a **3.2 mm** bayonet or trocar point pin and impact or screw into place. Remove the lateral drill bit and replace it with a pin in like fashion.


3) Now re-inspect the distal condyle contact locations of the ***PtoleKnee Surgical Guide***; the feet should still be in contact with the condyle surfaces. The anterior, medial, and lateral contact points above the epicondylar area should also be in contact. If needed, slight downward pressure on the anterior guide surface between the anterior pinholes may correct small contact gaps (less than 1.0 mm) (red circle) in the following picture.



4) Once all four contact point checks are completed, placement of the two anterior pins begins **(yellow circles)**. While still holding the Femur Guide in place, pre-drill the anterior lateral femur pin location, leave the lateral drill bit in place, and drill the medial side with a second drill bit. Replace the medial drill with a **3.2 mm** bayonet or trocar point pin and impact or screw into place. Remove the lateral drill bit and replace it with a pin in like fashion. In case a change of implant size is needed, use the different pairs of front pins accordingly as shown in the picture. Once all pins are inserted, the guide should be firmly locked in place. Examine the guide again once all pins are in place to verify that the guide has not skewed during pin placement.



5) With the guide fully pinned, and before cutting, examine the distal resection thickness on the medial and lateral sides using a resection checker (Angel-wing). The tool is standard in T.K.R. instrument sets. Check the potential cut carefully for thickness on both medial and lateral condyles for Varus/valgus angle and flex/extension angle. If potential alignment and cuts are as expected, proceed. The surgical plan provides the expected thickness of each condyle cut.

A kerf-less blade 110-120 mm x 1.27 mm is recommended.


Making the first cut:
Remove the lateral condyle distal pin **(red arrow)** and save it for reuse. The diagonal stabilizer **(blue arrow)** pin helps reduce possible guide shifting due to saw vibration. The medial pin **(black arrow)** remains in place until the lateral condyle cut is completed.

With the lateral femur condyle pin removed, resect the lateral condyle using an orthopedic saw.

If the area occupied by the diagonal stabilization pin requires resection, remove the pin and pass the saw blade across this area to complete the resection. When the lateral condyle is finished, remove the medial stabilization pin and reinsert carefully into the lateral location within the original hole to provide additional stability while cutting the medial condyle. Once both sides are cut, remove all pins and the guide. Check the flatness of the final cuts. Minor clean-up and smoothing of the cut surfaces may be undertaken if needed.


Sawing cortical bone can generate significant heat and may result in thermal necrosis. Irrigation of the saw blade with saline while cutting help reduce this heat. Multiple slow, smooth passes of the saw also help reduce heat and ensure a more accurate cut with less tendency to skive. When the femur cuts are completed, with the pins and guides removed, the resection appears similar to the photograph following.

Based on the implant system designated in the plan, the femur chamfer cuts may or may not follow now. Some systems immediately cut the tibia, followed by a resection gap check to assess gap dimensions and knee stability while corrections are still easily made.


Making the chamfer cuts:

Attach the chamfer-cutting (4:1) block after the revision of the femur is verified and found acceptable. Verification of this cut is crucial should conversion to standard manual instruments be necessary.

Correction of varus/valgus or flex-extension angles is difficult or impossible once the femur chamfer cuts are made. Occasionally, the distal holes can be challenging to identify following the initial distal cut. If this occurs, clean and irrigate the pinhole locations. Alternatively, before removing the femur guide, locate and re-clean the drill pinholes manually with a 3.2 mm drill bit by hand using the guide as an alignment tool.


Examine the Tibia Guide:

Shown are the 5 visible contact feet or points. Points 2, 3, 4 & 5 all contact the proximal surface of the tibial plateau. Point 1 contacts the anterior lip of the tibia. As with the femur guide, inspect all guide surfaces for damage, debris and that the guide serial number and patient identifying, and leg information is correct.


Preparation of the tibia:
Remove all meniscus tissue. Do not remove osteophytes from the tibia unless instructed to do so. Clean the bone around the A.C.L. stump and anterior to the medial tibial lip where the guide fits. Remove any soft tissue that may interfere with the proper placement of the ***PtoleKnee Surgical Guide***. Using a forked distractor/retractor posteriorly may be helpful to expose the tibia better and protect the Posterior Cruciate Ligament (P.C.L.) from additional injury when sawing

Placing the Tibia Guide:


Hold the tibia guide as shown above. The central pointer indicated (red arrow) aligns with the axis of the tibial spine when the guide is placed on the tibia. Orient the tibia guide with the central pointer aligned with the central spine of the tibia as shown (red arrow).


The ***PtoleKnee Surgical Guide*** contacts the anterior lip, and the cut slot extends medially on the tibia. It may be necessary to extend the original incision slightly to provide better guide access medially.


Place tibia guide against the anterior surface with the superior bilateral feet on top of the tibial plateau (pointer toward posterior) as shown. For reference, use the central arrow feature on the tibia guide to point along the central axis of the tibial spine, as shown above.


From the medial side, check the contact point of the medial guide feet within the circled area as above. Both feet should lightly touch the surface of the tibial surface along their length. The lateral side is checked in like fashion.


Also, check that the intended cut plane and natural tibial slope are approximately the same from the medial side, as shown above. The exception to this is when a tibial slope changes during planning to a different cut angle. In this case, verify that the new cut slope matches the planned cut slope.

While holding the tibia guide in place, install the two proximal tibia pins one at a time (collared and treated at the end), beginning with the lateral side (red arrow).

After the proximal pins are in place, continue to hold the tibia guide in place and re-verify all five (5) of the tibia guide contact points. If they are all still in contact with joint surfaces, the two 3.2 mm anterior guide pins are placed, beginning with the lateral pin the angled one from the top as first then the second one straight pin more medial.

Insert the anterior angled pin then insert the second anterior pin (more medial).
Remove one of the two tibia proximal pins (red circle) as the saw will work through the bone during the cutting phase.

Once the ***PtoleKnee Surgical Guide*** is securely anchored proceed with the alignment check.


Tibial alignment Check

A 1/4-inch drop rod from the standard instrument set may be placed as shown to verify it points to the center of the ankle.


Verification of the tibial cut-plain angle:

An estimate of the final Varus/valgus angle, posterior slope, and depth of cut of the tibial resection is conducted in the example shown below.


The tibial slope and depth of cut are expected to match the surgeons' plan.

The varus/valgus angle may be assessed using the long tibia axis and the guide cut plane. In this case, the plan called for zero degrees varus/valgus angle alignment.

If all verification checks match the plan or are as desired, the tibia may be cut with the power saw.

If the pinned guide does not replicate desired angles and slope, do not cut. Remove all the pins and tibia guide and continue with the remainder of the surgery using manual instruments.

Position the ***PtoleKnee Surgical Guide*** on the tibia as shown. Like the femur guide, careful checking of guide positioning before cutting is essential. Verify all contact anatomic contact points of the ***PtoleKnee Surgical Guide***. Verify guide stability. Check M/L thickness of the

proposed cuts. Check Varus/valgus angle matches the femur cut. Verify tibial slope matches the plan. If all is satisfactory, complete the tibial cut and proceed with the remainder of the T.K.R. procedure using the standard instruments.

Resuming the standard T.K.R. procedure:

For example, a typical next step verifies that the resection creates a rectangular resection gap between the distal femur and proximal tibia, which are flat and parallel in both flexion and extension.

Inserting a Flex/Extension gap checker from the implant tray allows assessment of the soft tissue tensioning and the final thickness of bone resected in both flexion and extension. Some implant brands assess fit during a "trial reduction." ***PtoleMedic System*** planning should result in a gap between the two cut surfaces allowing for the thickness of the metallic knee components and the thinnest polyethylene tibial insert. The minimum insert thickness varies by implant brand selected by the physician when setting their brand preference but is typically between 10 & 13 mm.

NOTE: *Some Implant systems and surgeons choose to perform an extension gap assessment first before completing the chamfer cutting of the femur. In this case, once the initial extension gap is verified as acceptable, the chamfer cuts are made, and the flexion gap is assessed.*

PtoleMedic System

Instructions for Use - Hospital Cleaning, Sterilization, and Disposal

Warnings:

- The PtoleMedic System has not been evaluated in a pediatric population therefor performance in such cases is unknown.
- The PtoleMedic System is a prescription only medical device.
- The PtoleMedic System is not a substitute for critical thinking and intra-operative adjustment of surgical goals based on education, training and experience of the practitioner.
- PtoleMedic System only provides, and documents useful alignment and orientation information based on specific individual anatomic data obtained from current MRI image sources.
- PtoleMedic System does not provide an absolute or only solution plan for joint replacement surgery, it only documents one possible approach, no surgical philosophy is recommended.
- These are patient specific, single use, disposable instruments.
- Do not attempt to reuse, recondition or re-sterilize.
- Do not alter the custom guides in any way.
- The surgical guides are to be used by a surgeon trained in the use of personalized instrument surgery (custom guides).
- Be aware that these patient specific instruments have been manufactured based of MRI scans of the patient. If the patient's anatomy has changed significantly since the time of the MRI scan, the instruments should not be used.
- The instruments should be properly cleaned and sterilized before use. Do not use if the instruments are broken, cracked or if debris is present.

- The surgical guides and their package are provided ***non-sterile only***.
- Federal Law (USA) restricts this device to sale by or on the order of a physician.

Cautions:

- The use of aged (>3 months) MRI image files is not recommended. Accuracy of planning and guide fit will diminish with evolving or changing disease processes.

Precautions:

- Use only MRI data of recent origin obtained per established PtoleMedic System designated MRI protocols.
- Care should be taken to minimize excessive heat buildup due to contact between the PtoleKnee Surgical Guides and other metallic instrumentation, such as drills. Excessive heat buildup can lead to deformation of the surgical guides.
- Do not place heavy instruments on top of the surgical guides during sterilization.

Limitations

- Metallic implants in or near the affected joint are known to interfere with the MRI images and may yield unreliable or useless images.
- The PtoleMedic System provides an estimate of implant sizing only. Exact implant size can only be determined during surgery and may differ from sizes projected during planning. Most estimated implant sizes will typically fall within one size of estimation.
- The PtoleMedic System is not for use in planning revision/replacement surgery in persons already having implants in the affected joints.
- Digital x-ray data is not acceptable for guide production the files must be MRI images.

Contra-indications

- Do not use in patients with active infection of the knee joint area.

NOTICE: Surgical Guides are intended to assist in the execution of the designated joint replacement surgery and can only be used in association with PtoleMedic System web software application. The guides are not reusable or transferable to any other person or surgery type.

Help Desk: If this document does not address your question, please contact Lento Medical Innovation, Inc. by:

Telephone: United States telephone support +1 (510) 413-3230
Internet: If you have access to the Internet, you can reach the Help Desk support team via:
World Wide Web at https://www.lentomedical.com/contact/

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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These services are available free of additional charges to all registered customers.

Indication for Use: The Surgical Guides are for clinical use as a template or guide for use in orthopedic surgery to assist the surgeon in selecting or positioning orthopedic implants and guiding the marking of tissue before cutting or pining for a specifically named patient.

CAUTION: United States Federal law restricts this device to sale by or on the order of a Physician

Acknowledgements:

Lento Medical Innovation, Inc. acknowledges the assistance of the following orthopedic surgeon for the expertise, guidance and time devoted to the development of the surgical protocol:

Benjamin Soo-Il Song, M.D., Diplomat American Board of Orthopaedic Surgery

Cleaning and Sterilization

The sonication cleaning process is more reproducible and reliable than manual cleaning methods. This method also affords the staff less exposure to the cleaning agents used. The cleaning method laboratory validation utilized a SonicWise Ultrasonics; Model: SW-308 sonicator.

Ultrasonic Cleaning

The preferred method of guide cleaning is ultrasonic (Sonicator) cleaning. This method is the most efficient and effective available today. Ultrasonic cleaning is up to 16 times more efficient than manual cleaning alone. Place instruments in the Sonicator for 15 minutes using a neutral pH enzymatic ultrasonic cleaning solution.

1) Before placing into the ultrasonic unit, rinse the instruments in warm running tap-water making sure to flush all channels and slots thoroughly.

2) Fill the Sonicator with fresh enzymatic cleaning solution* prepared according to manufacturer's instructions. A neutral pH ultrasonic solution should be used. Lento Medical Innovation, Inc. does not recommend a specific brand of enzymatic cleaner**.

3) Make sure instruments have plenty of room. Don't overload the ultrasonic cleaner. Sonicate the guides separately from all other instruments or articles.

4) Upon completion of the cycle, remove the guides immediately and rinse them under warm running tap-water as before.

5) Dry instruments thoroughly with a non-linting towel. Use compressed air to clear rinse water from all channels and slots. Ensure that no moisture remains on the Guides.

* ***Use distilled water or tap water with a neutral pH ultrasonic solution when filling the ultrasonic cleaner. Never use a manual soap in the ultrasonic cleaner and always follow proper dilution.***

** ***This cleaning procedure was validated using Miltex Surgical Instrument Cleaner (Phosphate free) neutral pH 7.05-7.65.".***

Cleaning validation was obtained using an FDA cleared cleaning agent.

<u>WARNING:</u> The Surgical Guides may not withstand exposure to automated cleaners operating at 285°F (141°C) and above or which use live-steam jets as cleaning features.

Use only FDA cleared sterilization pouches or wraps. Place each guide in a 7"x12" peel pouch - Mylar type (suggested)

Validated sterilizer parameters are:
Sterilizer Type: Steam, Vacuum (Pre-Vac) Sterilizer
Exposure (Holding time): 3 minutes
Temperature: 275°F or 135°C (EU)
Drying Time: 32 Minutes

After sterilizing, the device(s) should be allowed to cool to room/ambient temperature or until adequately cooled for safe handling.

Disposal

The Surgical Guides are constructed from recognized biocompatible materials and do not represent hazardous substances in terms of disposal. Once used, the devices must be considered biohazard due to their exposure to biological fluids of a patient. Disposal must be conducted in accordance with state local and hospital policy for handling biohazard materials.

In the event the devices are to be returned to Lento Medical Innovation and the devices are unused, no precautions are necessary and standard shipping may be used. However, if the devices have been used, they are considered biohazard and must be cleaned and sterilized before return. Refer to the "Cleaning and Sterilization" section of this document for instruction on how to clean & sterilize the device for return.

PtoleMedic System

Instructions for Use – Italian Translation

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