PtoleMedic System

Software Instructions for Use

Document ID: LB-72-01-007
Model: Revision H

October 2, 2026

Author:

Name: **Alex Park** **10/2/2026**
Date

Reviewed and Approved by:

Name: **Won Choi** **10/2/2026**
Date

TABLE OF CONTENTS

1. Instructions for Use - Adding a New Case…… 7
2. Instructions for Use - Uploading MRI to Database (MRI Center)……15
3. Instructions for Use - Template Design Approval (Surgeon)……21
4. Instructions for Use – Clinical MRI Scanning……31
5. Instructions for Use – MRI Protocol Quick Reference Guide……48

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Introduction

PtoleMedic System *web-based software* offers a simplified Total Knee Replacement (TKR) surgery planning interface. Surgical planning software is widely available but is generally PC or local server based. The use of such programs tends to be time consuming, expensive and requires frequent updating and/or IT maintenance. The PtoleMedic System offers **"log-on-and-use"** surgical planning tools and custom made cutting guide creation. And the system is **low-to-no maintenance** requiring no IT support to install, use or maintain. The PtoleMedic® System web-based software provides direct, focused assistance on setting angles, rotation, cut plane, offsets and implant sizing as appropriate to TKR surgery.

Warning:

- The PtoleMedic System has not been evaluated in a pediatric population therefore performance in such cases is unknown and not recommended.
- The PtoleMedic System software is a prescription only medical device.
- The software is not a substitute for critical thinking and intra-operative adjustment of surgical goals based on education, training and experience of the practitioner.
- Surgical planning software only provides and documents useful alignment and orientation information based on specific individual anatomic data obtained from current MRI image sources.
- The software does not provide an absolute or only solution plan for joint replacement surgery, it only documents one possible approach, no surgical philosophy is recommended.
- The use of the software on mobile platforms, such as Google Android, Apple iPhone, or similar tablet-type devices has not been validated and is not recommended at any time.
- These are patient specific, single use, disposable instruments.
- Do not attempt to reuse, recondition or re-sterilize.
- Do not alter the instruments in any way.
- PtoleMedic System guides are to be used by a surgeon trained in the use of personalized instrument surgery (custom guides).
- Be aware that these patient specific instruments have been manufactured based of MRI scans of the patient. If the patient's anatomy has changed significantly since the time of the MRI scan, the instruments should not be used.
- The instruments should be properly cleaned before sterilization. Do not use if the instruments are broken, cracked, or when loose powder is present.
- The instruments in the package are provided non-sterile.
- The instruments in the package must be sterilized prior to use.
- Federal Law (USA) restricts this device to sale by or on the order of a physician.

Cautions:

- The use of aged (>3 months) MRI image files is not recommended. Accuracy of planning and guide fit will diminish with evolving or changing disease processes.

Precautions:

- Use only MRI data of recent origin obtained per established PtoleMedic System designated MRI protocols.

- Care should be taken to minimize excessive heat buildup due to friction between the PSI instruments and other instrumentation, such as drills. Excessive heat buildup can lead to deformation of the PSI instrument.
- Do not place heavy instruments on top of the cutting guides instruments during sterilization.

Limitations

- Metallic implants in or near the affected joint are known to interfere with the MRI images.
- Provides estimate of implant sizing only. Optimal implant size can only be determined during surgery and may differ from sizes projected during planning. Most estimated implant sizes will typically fall within one size of estimation.
- The software is not for use in planning revision/replacement surgery in persons already having implants in the affected joints.
- DICOM to JPEG file creation – The use of planning images based on the JPEG image format limits some file attributes encoded with the original DICOM images. For instance, some information contained in the DICOM header including image manipulation functions like pan, zoom, full gray scale control, contrast adjustment and pixel resizing are limited or unavailable.
- Digital x-ray data is not acceptable for planning purposes; the files must be MRI images.

Contra-indications

- Do not use in patients with active infection of the knee joint area.

User System Requirements

1. High speed internet access
2. A color monitor or display with VGA image display capabilities with 1280 x 800 resolution (recommended)
3. A mouse, touch pad or similar input device to interact with the surgery case
4. A keyboard (virtual or physical) for text input
5. A PC/Workstation w/1GB of RAM or more capable of running Windows or Mac OS
6. A Web Browser, (Internet Explorer, Firefox, Google Chrome or Safari)

Clinical Benefit

PtoleMedic System enhances the accuracy of implant placement, potentially leading to better outcomes in surgery and improvement in quality of patient life.

Indication for Use
The PtoleMedic System is an on-line orthopedic surgical planning software system. MRI images supply data sufficient to allow accurate modeling of anatomy for on-line surgical planning before knee joint replacement surgery. The surgeon preoperatively plans, reviews, adjusts and orients implant images to select implant size and create an idealized surgical plan for the first bony cuts only.

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively

Intended Users:
These instructions are intended to guide Surgeons, Physicians, and MRI Technologist in using the PtoleMedic System software for creating new patients, MRI scanning of patients, Uploading patient MRI scans, and approving template designs.

NOTICE: ***Software Anomalies***. The PtoleMedic System software supports the following Web Browsers, Internet Explorer (Ver 9.0) or newer, Firefox (Ver 15.0) or newer, Google Chrome (Ver 25.0) and Safari (Ver 4.0) or newer. Such browsers may be used with various operating systems, including those supporting Windows OS (XP, 7, Vista) or newer or Mac OS X (Ver 10.0) or newer.

Depending on your browser or operating system your screen may look slightly different depending on the software environment. This may include available fonts, spacing between graphical user interface (GUI) components, and some other minor appearance items.

NOTICE: The PtoleMedic System Web software does not interact directly with DICOM generating equipment (PACS devices). The image files needed by the PtoleMedic System Web software originate from the MRI/CT center as DICOM files which are transferred to the company server by Imaging Center personnel using DICOM compliant processes. File transfer protocols or file transfer requests are not generated or provided.

NOTICE: The PtoleMedic System Surgical Planning Software program does not interact with or affect DICOM image source equipment (PACS devices).

NOTICE: In the event of a serious incident in relation to the PtoleMedic System, report the serious incident to Lento Medical Innovation, Inc. and the competent authority of the Member State in which the user and/or patient is established.

Updates: For updates to this document, please contact the Help Desk support team.
Help Desk: If this document does not address your question please contact Lento Medical Innovation, Inc.

Telephone: The Help Desk and Customer Service support teams can be reached via:
United States and Canada during normal operating hours +1 (510) 413-3230.

Company URL: http://www.lentomedical.com

PtoleMedic System URL: http://www.lentomedical.net

MD

UDI

Basic UDI-DI: ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that the use of MRI scanners in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be performed, please verify the MRI order with the physician.

WARNING: Making custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. For this reason, it is generally recommended that such patient not be entered into the patient database. Please verify this request with the ordering physician.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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PtoleMedic System

Instructions for Use - Adding a New Case

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively

Intended Users:
This is a reference guide for the office manager, surgical planner/scheduler or physician as an aid in smoothly registering a new patient into the PtoleMedic System database.

MD

UDI

Basic UDI-DI: ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that the use of MRI scanners in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be performed, please verify the MRI order with the physician.

WARNING: Making custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. For this reason, it is generally recommended that such patient not be entered into the patient database. Please verify this request with the ordering physician.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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

Introduction

These instructions are intended to guide the surgeon or surgeons' representative in creating a new case for a patient. Adding a new case through the Lento Medical Innovation web portal (www.lentomedical.net) will create a new order for the cutting guide.

1. **Log In**
 Lento Medical Admin or Sales personnel will assist the surgeon with registering an account for the Lento Medical Innovation web portal. Once a username and password are created, the surgeon will use the information to log in.
 i. Enter Username
 ii. Enter Password
 iii. Click "Login" button


2. **Surgeon Landing Page**
 The surgeon landing page will appear once logged in.

3. Adding a New Case

From the Surgeon Landing page, click the “NEW ORDER” button to create a new case.


As seen in the above image, select surgery type by clicking “TKR”.

If the Surgeon User already has more than one pre-set preference, then user can select a preference (maximum of 3 preferences).

4. Entering Information for New Case

A. Add New Patient

A new page will be presented to add the following information:

I. First Name*
II. Middle Name
III. Last Name*

IV. Sex*
V. Date of Birth*
VI. Date of Surgery*
VII. Surgery Side*


Note: * fields are mandatory

B. Check Preference Information

I. Check distributor information
II. Check Hospital information
III. Check MRI Center information
IV. If any of the information are incorrect, STOP contact your distributor or Lento admin.
V. If you want to change the preference, go back to the previous screen and select a different preference.
VI. If all the information is correct, move forward to the next step.

C. Confirm Surgery Preference Information

I. Billing

i. Billing party pre-set, (Hospital, Distributor, or Surgeon)

II. Shipping
 i. Ship to location preset (Hospital, Distributor, or Surgeon)

III. Implant
 i. Implant Brand pre-set.

IV. Alignment Default Position
 i. Select alignment default position (Mechanical or Natural Boundary).

V. Click “Submit”, to submit new order for case creation.


If you want to change the preference, go back to the previous screen and select a different preference.

5. New Case Creation

Once a case is created, this will create a new order in the PtoleMedic System based on the information provided. A unique serial number, the “case ID number”, is assigned to the new patient case as seen below.


Sample Case:

If the new case created requires any changes or updates, please contact a Lento Medical personnel or the sales representative for further assistance.



I. Collapse button will hide preference information to minimize window space.

6. Tutorial Video
 i. You can also access a tutorial video link located at the top right menu for a visual demonstration of how to create a new order.


7. Logout

If the surgeon user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.

PtoleMedic System

Instructions for Use - Uploading MRI to Database (MRI Center)

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:
This is the reference guide for the MRI Technologist as an aid in uploading the MRI scanning images to the Lento Medical Innovation cloud database. There is a separate full MRI reference guide for use in positioning and patient orientation. Also, a quick MRI setting reference guide is available for use by experienced technologists and both are available upon request if they have not already been supplied or is misplaced.

MD

UDI

Basic UDI-DI: ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery be attempted.

WARNING: Federal Law (USA) restricts this device to sale by or on the order of a physician.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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

Introduction

These instructions are intended to guide the MRI Center to upload MRI files through the Lento Medical Innovation web portal (www.lentomedical.net).


1. **Log In**
 Lento Medical Admin or Sales personnel will assist the MRI representative with registering an account for the Lento Medical Innovation web portal. Once a username and password are created, the surgeon will use the information to log in.
 i. Enter Username
 ii. Enter Password
 iii. Click “Login” button

2. **MRI Center Landing Page**
 User landing page will appear listing all surgeons. See below.


Click “Surgeon Orders” button to see a list of surgeons that have cases waiting for MRI files to be uploaded, as seen below. User will only see Surgeons in the list that have open orders waiting for MRI Files.


3. MRI Uploading Page

From the Surgeon list, User will the click the surgeon's last name, which will direct the User to Case List page that shows all cases that are waiting for MRI files to be uploaded for the respective surgeon. The status of the case will show as "Ready to Upload MRI", see below image. NOTE: Before uploading MRI files, verify the patient information and the Case ID (unique serial number):

User will click the case number and the MRI upload page will appear, see below.

Note: before uploading the MRI file, put all the MRI files "Coronal Kee", "Coronal Hip", and "Coronal Ankle" into one folder.

4. Uploading MRI

There are two ways you can upload MRI files.

A. Drag and drop
 User can drag the MRI file folder and drop the file folder into the "Drop Zone"

B. Manual folder selection
 User can click the "Drop Zone" to select the folder and click the "Upload" button.


5. **Uploading Progress**

Once the User as begun uploading the MRI files, there will be an uploading progress bar as seen below.

6. **MRI upload confirmation**

Once MRI files are uploaded correctly and successfully, a confirmation popup will present itself as seen below.


Once you click “OK”, user will be re-directed to user landing page.

Note: once MRI is uploaded successfully, case will no longer be on the case list.

Repeat from Step 2 to upload MRI files for a different case.

7. **Logout**

If the user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.



8. Follow up

Our technican will review the MRI files that had being upload and contact you if there is any issues.

9. Tutorial Video

i. You can also access a tutorial video link located at the top right menu for a visual demonstration of how to upload MRI Files.

PtoleMedic System

Instructions for Use - Template Design Approval (Surgeon)

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:
These instructions are intended to guide the surgeon in the surgical planning and template design, review and approval stage for an order of the PtoleMedic System guides. Approval is provided online through the Lento Medical Innovation web portal (www.lentomedical.net).

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

Caution: Please do not share your Lento login information and keep your password confidential. Please reset your password or contact Lento Medical Administrator (support@lentomedical.com) if you have trouble logging into the system.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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
1. Log In

Lento Medical Admin or Sales personnel will assist the surgeon with registering an account for the Lento Medical Innovation web portal. Once a username and password is created, the surgeon will use the information to log in.

i. Enter Username
ii. Enter Password
iii. Click “Login” button


2. Surgeon Landing Page

The surgeon landing page will appear listing basic information provided by the Surgeon. To see active cases, press the “Current Orders” button.


The status of a case that requires surgeon approval will show “Surgeon Approval Ready”


3. View Surgeon Planning Reference

The Surgeon can click the case ID number (unique serial number) and the patient case detail page will appear.

From the patient case detail page, surgeon can view case information.

Press “EXPNAD” button to view preference details shown below.


Press “COLLAPSE” to hide the preference.


Press the “Surgeon Approve” button and the Surgeon will be directed to the “Surgeon Template Design Approval” Page

4. Surgeon Template Design Approval Page

On this page, surgeon will have the option to view each coronal, sagittal, and axial view of the MRI surgical planning results.


Each Coronal, Axial, and Sagittal view will be enlarged by double clicking the corresponding MRI image, shown as below.


Each arrow button will move the implant positioning (left, right, up, and down)
The rotate icon/button will rotate the implant by 3 degrees.


The overall limb alignment may be adjusted by clicking the right or left arrow on either side of the knee graphic image. The circle at the top of the graphic represents the center of the femoral head and the brackets at the bottom represent the center of the ankle joint. Each right or left arrow click adjusts the varus/valgus angle of the knee 0.5 degree or 0.5mm, max +- 5 degree or mm. The knee image is moved medially or laterally while maintaining the upper and lower boundary area within the femoral head and the median of the ankle. The resection lines of the rectangular resection box remain parallel to the floor.


Hip button will display Femoral Head Center in pop up.

Femur button will display Distal Femur Center in pop up.

Tibia button will display Proximal Tibia Center in pop up.

Ankle button will display Ankle Center in pop up.


After the surgeon has finished reviewing the planning data, at the bottom there will be four action buttons displayed. A button labeled reset, a button labeled reject, button labeled accept, and a button labeled close.

If the user presses the close button, system shall discard any saved changes and restore the original values and close the surgeon template approval page. The page will be redirected to the case detail page.

If the user presses the reset button, all the values will revert back to the original values from the beginning prior to any adjustments made by the surgeon user.

If the user presses the reject button, the system shall reject the entire plan. The surgeon will need to submit a reason for rejecting the plan as shown below.


Press Ok to submit the reason for rejection.


If the user presses the accept button, the OMM shall accept the current planning values. The user shall not be allowed to change the values after acceptance.

5. Tutorial Video

i. You can also access a tutorial video link located at the top right menu for a visual demonstration on how to approve the Template Design.


6. Logout

If the surgeon user is finished with their session, they can logout by going to the top menu and click logout. The user will be redirected to the login page.

PtoleMedic System

Instructions for Use – Clinical MRI Scanning

This is a quick reference guide for the MRI Technologist as an aid in rapidly setting up the listed MRI scanning equipment. There is a full MRI reference guide for use in positioning and patient orientation which is available upon request if it has not already been supplied or is misplaced. There is also an MRI upload reference guide to assist with uploading completed images to the Lento Medical Innovation cloud database.

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery not be attempted.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer of Record:
Kalms Consutling GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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A. PATIENT SET-UP

1. Patient should be metal/implant free in the scan volume.
2. Place patient supine, feet first, for lower extremity with affected limb as close to IsoCenter as possible
3. Consider using small foam wedges to stabilize or support the joint.
4. Use tape, sponges and/or Velcro straps on extremity to maintain anatomic “toes up” position.
5. Remind patient to remain motionless during the entire scan.
6. Landmark or Zero table
7. Begin scanning at the Knee series and move to Distal Tibia (Ankle) for the Coronal T1 Series and complete the study with the Hip series (see Figure 11).
8. Phased array coils (receive only coils) or the “Body Coil” may be used to scan the hip and ankle without removing the knee surface coil. New GEM suite GE systems have built in phased array coils in table.


B. SCANNER SET-UP

1. Enter patient name (Last, Middle, First)
2. Enter patient DOB. (dd/mm/yyyy)
3. Enter scan date (dd/mm/yyyy)
4. Enter left or right joint (Right/Left)
5. Enter patient gender (M/F)
6. Enter name of Surgeon (Last, Middle, First)
7. Enter name of imaging center

C. CORONAL FSE PD KNEE MRI PROTOCOL

1. Perform one 3 planes (Axial, Sagittal, Coronal) locator series.
2. Refer to "PtoleKnee Imaging Reference for 1.5T and 3T Scanners section for the Coronal scan parameters

D. CORONAL FSE PD KNEE SLICE POSITIONING

(a) **(b)**

Figure 2
Slice Positioning through Posterior Reference and Distal Reference Lines

1. Using the best axial and coronal locator images that show in Figure 2
 i. both posterior lateral and medial condyles, place slice parallel through the posterior reference for correct slice angle
 ii. Both distal lateral and medial condyles, place slice parallel through the distal reference for correct slice angle.


Figure 3
Slice Coverage

2. Slice coverage must include all the Femur including the condyles, Tibia and Patella as shown in Figure 3.
3. Adjust T1 T2 setting to show optimal image quality

Figure 4
Centering FOV

4. Using the best sagittal locator image that approximately shows the knee center, adjust the FOV to be centered at the “point” at the bottom of the distal femur as shown above in Figure 4.




Figure 5
Axial axis of the Knee joint

Not Acceptable

Acceptable

5. Find the true Axial plane or axis of the knee joint as shown in Figure 5 by the yellow dash line. In addition, locate the sagittal MRI scout scan slice representing approximately the mid tibia, and tilt about 30 degrees. The upper bound of femur coverage must include femur anterior cortex or the beginning of the femoral shaft. The lower bound of tibial coverage must include the tibial shaft as much as possible, shown in Figure 6.


Figure 6
Slice Coverage

6. **Enter File Series Name: “Coronal Knee.”**

A. General Electric

3plane Localizer Scan Parameter	GP FLEX (GE Users), 4mm x 1mm Skip 24cm FOV, Matrix 256x192
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL
Mode	2D
Imaging Option	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)
TE (Echo Time)	Min Full
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	18cm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

B. Philips

3plane Localizer Scan Parameter	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x 192
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	Fold Over, Suppression ON, 100% Sampling, "Default" Selected for Distortion Correction
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	8
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

C. Siemens

3plane Localizer Scan Parameter	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	100% Phase Oversampling 3D Distortion Correction Filter ON
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700**)
Flip Angle (deg)	120
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	180mm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor	0%
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	H/F

***Higher resolution than 256 x 256 is allowed, using software.**
256 x 256 is optimal and preferred.

D. Toshiba

3plane Localizer Scan Parameter	4mm x 1mm Spacing, 24cm FOV, Matrix 256 x 192
Coronal Knee: Pulse Sequence	TSE CORONAL
Mode	2D
Imaging Option	Fold Over Suppression ON, IDC (Intelligent Distortion Correction) Selected
TE (Echo Time)	~24 to 35 (28 nominal)
TR (Repetition Time)	**1600 average (min 1500 to max 1700)**
Flip Angle (deg)	90, Flop Angle 160
Echo Train Length (ETL) **Turbo Spin Factor** (TSF)	7
FOV (Field of View)	18cm
Slice Thickness(mm)	2.5mm
Spacing/Skip/Gap/Distance Factor(mm)	0mm
Scan Matrix/Voxel Size (Base Resolution x Phase Resolution)	256 x 256
NEX/NSA/Averages	2
Frequency Direction	S/I

***Higher resolution than 256 x 256 is allowed, using software.**
256x256 is optimal and preferred.

E. CORONAL T1 DISTAL/ANKLE MRI PROTOCOL

1. Perform one 3 plane (Axial, Coronal, Sagittal) locator scan
2. Perform a Coronal (T1FSE) series using the minimum scan
 i. parameters:

Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	10
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Ankle"


Figure 1
Coronal Slice Placement

F. CORONAL T1 DISTAL TIBIA/ANKLE SLICE POSITIONING

View the **Sagittal** locator to place the slice Coronally through the long axis of the Tibia. *The Axial locator slices are not used for alignment.* Center the FOV approximately one inch above the Tibial Talar joint space. See Figure 1 for slice position and FOV centering.

G. CORONAL T1 HIP MRI PROTOCOL AND SLICE POSITIONING

1. Perform one 3 plane (Axial, Sagittal, Coronal) locator scan.
2. Perform a Coronal (T1 FSE) series using the minimum scan parameters total below:
3. Place Coronal slice plane parallel through long axis of the femoral neck when viewing the axial locator as shown in Figure 8. Center the 14 slices to cover the Femoral Head, Neck, and Greater Trochanter as shown in Figure 9.


Figure 10
Center Hip FOV

Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	14
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Hip"

I. QUALITY CHECK AND METAL ARTIFACT INSTRUCTIONS

a. Coverage: Include 1 or 2 slices of soft tissue on both medial and lateral side of bony knee anatomy.
b. Motion: Make sure sequential transition of the bone/cartilage outline does not jump or shift back and forth between slices
c. If patient presents with metal in the affected lower extremity, follow these instructions
 - Metal Hardware Knee – Do **NOT** Scan the patient. Refer them back to the surgeon
 - Metal Hardware Hip and/or Distal Tibia – Use MARS (Metal Artifact Reduction Sequence) series. Keep in mind the FOV, Slice Thickness, and Spacing/Gap/Skip/Distance Factor must be set to our parameters. Do not use FAT SAT at any time.
 - For those facilities that do not use a MARS series, increase the bandwidth and NEX/NAQ/NSA to help increase signal to noise ratio.

II. FILE FORMAT & OUTPUT

a. Save all image series in DICOM format
b. Upload the following DICOM files to Lento Medical Innovation's database (www.lentomedical.net)
 i. **Coronal Knee**
 ii. **Coronal Hip**
 iii. **Coronal Ankle**




Figure 11
Scanning Order

I. KNEE, ANKLE, AND HIP MRI SAMPLE IMAGES


a. KNEE MRI – file name: Coronal Knee


(a)

(b)

(c)

(d)

(e)


(f)


(g)

b. **ANKLE MRI – file name: Coronal Ankle**


(a)

(b)


(c)

(d)


c. **HIP MRI – file name: Coronal Hip**


(a)

(b)

(c)

(d)

PtoleMedic System

Instructions for Use – MRI Protocol Quick Reference Guide

Lento Medical Innovation, Inc.
Imaging Quick Reference Guide
for 1.5T and 3T Scanners

Intended Purpose:
The PtoleMedic System is intended to be used as an orthopedic surgery planning software and patient specific surgical instrument design and production software to assist the surgeon with the sizing and anatomic positioning of Total Knee Replacement implants (TKR) components intraoperatively.

Intended Users:
This is a quick reference guide for the MRI Technologist as an aid in rapidly setting up the listed MRI scanning equipment. There is a full MRI reference guide for use in positioning and patient orientation which is available upon request if it has not already been supplied or is misplaced.

MD

UDI **Basic UDI-DI:** ++G283PtoleMedicSystemFC
UDI-DI/PI: +G283PMS0010/&&71E
SRN: US-MF-000028244

WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.

WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery be attempted.

US

Manufacturer:
Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230

Importer:
Kalms Consulting GmbH
Rheinstr. 45-46
12161 Berlin, Germany

The following organization is the Authorized Representative for the PtoleMedic System:

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3T MRI Scanners	General Electric	Philips	Siemens	Toshiba
3 plane Localizer Scan Parameters:	**GP FLEX (GE Users)** 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution	4mm x 1mm Spacing 24cm FOV, Matrix 256x192
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL	TSF CORONAL
Mode	2D	2D	2D	2D
Imaging Options	No Phase Wrap ON, Tailored RF, 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression ON, 100% Sampling, “Default” Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON	Fold Over Suppression ON, IDC (Intelligent Distortion Correction) Selected
TE (Echo Time)	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
TR (Repetition Time	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**
Flip Angle (Deg)	90	90	120	90 Flop Angle 160
Echo Train Length (ETL), Turbo Spin Factor (TSF)	7	8	7	7
FOV	18cm	180mm	180mm	18cm
Slice Thickness (mm)	2.5mm	2.5mm	2.5mm	2.5mm
Spacing/Skip/Gap/Distance Factor (mm)	0	0	0%	0
Scan Matrix/Voxel Size Base Resolution x Phase Resolution	256 x 256	256 x 256	256 x 256	256 x 256
NEX/NSA/Averages	2	2	2	2
Frequency Direction	S/I	S/I	H/F	S/I

3T MRI Scanners	General Electric	Philips	Siemens
3 plane Localizer Scan Parameters:	**GP FLEX (GE Users)** 4mm x 1mm Skip, 24cm FOV, Matrix 256 x192,	4mm x 1mm Gap, 240mm FOV, Voxel Size to equal 256 x192	4mm x 25% Distance Factor, 240mm FOV, Base Resolution 256 x 80% Phase Resolution
Coronal Knee: Pulse Sequence	FRFSE-XL CORONAL	TSE CORONAL	TSF CORONAL
Mode	2D	2D	2D
Imaging Options	No Phase Wrap ON, TRF (Tailored Radio Frequency), 3DGR (3D Geometry Correction for software version 23+)	Fold Over Suppression R/L, 100% Sampling, "Default" Selected for Distortion Correction	100% Phase Oversampling, 3D Distortion Correction Filter ON
TE (Echo Time)	Min Full	~24 to 35 (28 nominal)	~24 to 35 (28 nominal)
TR (Repetition Time	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**	**1600 average (min 1500 to max 1700)**
Flip Angle (Deg)	90	90	120
Echo Train Length (ETL), Turbo Spin Factor (TSF)	7	8	7
FOV	16cm	160mm	160mm
Slice Thickness (mm)	2.5mm	2.5mm	2.5mm
Spacing/Skip/Gap/Distance Factor (mm)	0	0	0%
Scan Matrix/Voxel Size Base Resolution x Phase Resolution	256 x 256	256 x 256	256 x 256
NEX/NSA/Averages	2	2	2
Frequency Direction	S/I	S/I	H/F

PtoleMedic System

Instructions for Use – Italian Translation

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MDSS GmbH